
Produktinformation

ELEC[®]PLUS

Kugelpf-Hüftgelenkimplantate
Femoral head implants

ARTIQO

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Anatomy as a benchmark

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We develop, produce and market artificial joints based on the needs and scientific findings of trauma surgeons and orthopedists. Close networking with scientists and clinically active surgeons is focused on the well-being of patients. The already good clinical results of joint replacement are to be further improved through consistent further development in the areas of materials research, implementation of new mechanical and biomechanical findings, improvement and refinement of surgical techniques, education and training of hospital personnel, and perfection of logistics.

ELEC®plus CERAMIC FEMORAL HEAD

Product features

- | the plus in safety due to 50% higher fracture resistance compared to conventional aluminum oxide heads.
- | extremely abrasion-resistant, bioinert and mechanically highly resilient
- | mixed-phase, ceramic high-performance material (zirconium oxide-reinforced aluminum oxide)
- | available as standard in three plain bearing sizes (28, 32, 36) and three neck lengths (S, M, L) and optionally in XL for sizes 32 and 36
- | also available as revision option with metal sleeve
- | no chromium oxide in the material
- | meets the requirements of the ISO 6474-2 standard

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Femoral and acetabular implants for hip arthroplasty made of ELEC®plus high-performance ceramic materials have proven their worth in more than 10 years of clinical use.

As a low-abrasion alternative to metal femoral heads in combination with polyethylene articulation inserts, ELEC®plus femoral heads can be combined with all common stem materials. ELEC®plus femoral heads contribute to the safe fitting of modern endoprosthesis systems due to their high fracture strength and extremely wear-resistant material. For revision cases, we offer ELEC®plus revision femoral heads that can be mounted on a stem that remains in situ.

ELEC®plus components are made from the purest starting materials without the use of coloring admixtures (e.g. chromium) and are completely harmless with regard to remaining in the human body for decades.

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

Controlled production processes are the basic prerequisite for the manufacture of ceramic implants according to specifications.

Based on more than three decades of experience in powder preparation, thermal process technology and finish machining, this requirement is met in production.

State-of-the-art manufacturing processes guarantee reliable production of the ELEC®plus implants. A high degree of automation, customized process solutions and suitable clean room environments ensure reproducible component quality and thus patient safety.

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
Continuous monitoring of all manufacturing steps as well as comprehensive approval testing of every single component assure a maximum of product safety.

The outstanding biocompatibility of and ELEC®plus implants has been – beyond obligatory testing – validated in long-term human cell and tissue culture and is completely confirmed in ten years of clinical use.

All processes are controlled and monitored via quality management in accordance with ISO 13485.

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Abrasion resistance: Simulator test according to ISO 14242-1 and ISO 14242-2

Articulation	Components	Abrasion femoral head	Acceptance criterion according to ISO 6474-2 and comparison articulation
Ceramic/ceramic (aged according to ISO 6474-2)	ELEC®plus femoral head / ELEC®plus cup insert	0,008 mm³/ Mio	5mm³/Mio
Ceramic/PE (aged according to ISO 6474-2)	ELEC®plus femoral head / crosslinked polyethylene cup insert	0,011 mm³/ Mio	5mm³/Mio

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Burst load and cyclic fatigue in accordance with ISO 6474-2 and ISO 7206-10

Components	ELEC®plus femoral head	Determined forces (MW, n=5)	Acceptance criterion (MW, n=5)
Burst load worst case Femoral head component	ELEC®plus 28L femoral head on CoCr stem taper 12/14	69 kN, min: 59 kN (MW, n = 5)	Mean Value ≥ 46 kN, Min-Value: 20 kN, out of 5 samples **
	ELEC®plus 28L femoral head on Ti-stem cone 12/14	85 kN, min: 59 kN (MW, n = 18) *	Mean Value ≥ 46 kN, Min-Value: 20 kN, out of 5 samples **
Burst load worst case Femoral head component post fatigue	ELEC®plus femoral head / on CoCr stem taper (ISO 5832-4)	53 kN (n=3; Min-Value:43 kN) *	Min-Value: 20 kN after 10 ⁷ cycles (1.4 to 14 kN) of dynamic load **

Abzugtest – ASTM F2009

Pull-off test - ASTM F2009

* Values for specific combination of stem material, stem taper geometry and taper surface roughness, other stem materials, stem taper geometries and taper surface roughness might result in higher or lower values in the described tests

Components	ELEC®plus femoral head	Determined forces (MW, n=5)	Acceptance criterion (MW, n=5)
Pull-off force worst case connection	ELEC®plus 28L femoral head on CoCr stem taper 12/14	1.459 *	>1.000 N

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Material composition

* Values for specific combination of stem material, stem taper geometry and taper surface roughness, other stem materials, stem taper geometries and taper surface roughness might result in higher or lower values in the described tests

Property	unit	Specification ISO 6474-2 Type X	Specification ELEC®plus
Average relative bulk density	%	≥ 99	≥ 99
Chemical composition			
Alumina Matrix, Al ₂ O ₃	Ma.-%	60 - 90	75*** - 80
Zirconia ZrO ₂ + HfO ₂	Ma.-%	10 - 30	25*** - 30
Amounts of HfO ₂ in ZrO ₂	Ma.-%	≤ 5	≤ 3
Intended additives	Ma.-%	≤ 10	≤ 3
Limits of impurities	Ma.-%	≤ 0,2	≤ 0,2

** FDA Guidance Document for the Preparation of Premarket Notifications for Ceramic Ball Hip Systems, January 1995

*** Nominal value: small deviations permitted