ELEC[®]plus Revision

Purity + Strength = Permanence

Instruction for use

HiPer Medical AG
 Ziegeleistrasse 7
 16727 Oberkrämer OT Eichstätt
 Germany
 hiper@hiper.de
 www.hiper-medical.de
 Tel.: +49 (0) 3304 3977-0
 Fax.: +49 (0) 3304 3977-90

DE	Gebrauchsanweisung ELEC®plus-Revisionshüftkopf mit Titanhülse
BE	Рэвізійныя сцегнавыя галоўкі ELEC®plus з тытанавай гільзай
	Инструкция за използване
BG	Ревизионни бедрени глави ELEC®plus с титанова втулка
	Инструкцыя па выкарыстанні
CS	Návod k použití ELEC®plus – revizní hlavice stehenní kosti s titanovým pouzdem
DA	Brugsanvisning ELEC®plus-revisionshoftehoveder med titanmuffe
EL	Οδηγίες χρήσης Κεφαλές ανάταξης άρθρωσης ισχίου ELEC®plus με χιτώνιο τιτανίου
EN	Instructions for use ELEC®plus revision femoral heads with titanium sleeve
ES	Instrucciones de uso Cabezas femorales de revisión ELEC®plus con casquillo de titanio
ET	Kasutusjuhendit Revisjoni reieluupead ELEC®plus koos titaanhülsiga
FI	Käyttöohjetta ELEC®plus-revisioreisiluunpäät titaaniholkilla
FR	Mode d'emploi Têtes fémorales de reprise ELEC®plus avec manchon en titane
HR	Upute za korištenje Revizijske femoralne glave ELEC®plus s titanijskom čahurom
HU	Használati utasítás ELEC®plus-revíziós combcsontfejek titánhüvellyel
IT	Manuale d'istruzioni Teste femorali da revisione ELEC®plus con manicotto in titanio
LT	Naudojimo instrukcijos „ELEC®plus“ kontrolinės šlaunikaulio galvutės su titano įvove
LV	Lietošanas instrukcija ELEC®plus revīzijas gūžas galviņas ar titāna uzstavu
NL	Gebruiksaanwijzing ELEC®plus-revisieheupkoppen met titanium huls
NO	Bruksanvisning ELEC®plus-revisjonshoftehode med titanhylse
PL	Instrukcja użycia Głównki rewizyjnej ELEC®plus z tuleją tytanową do endoprotezy biodra
PT	Instruções de utilização Cabeças do fémur de revisão ELEC®plus com manga de titânio
RO	Instructiuni de utilizare Capete femurale de revizie ELEC®plus cu manşon din titan
RU	Инструкция по применению ревизионных бедренных головок с титановой втулкой ELEC®plus
SK	Návod na použitie Revízne hlavice bedrového kĺbu ELEC®plus s titánovou objímkou
SI	Navodila za uporabo Revizijske glave stegenice s tulcem iz titana ELEC®plus
SV	Bruksanvisning ELEC®plus revisionshöft huvuden med titanhylsa
TR	Kullanım Talimatı Titanyum kovanlı ELEC®plus revizyon femur başları
UK	Голівка ревізійного ендопротеза ELEC®plus з титановою гільзою
	Інструкція із застосування
VI	Hướng dẫn sử dụng Đầu xương đùi loại Revision ELEC®plus với lớp bọc titan
ZH	使用说明 带钛合金套筒的 ELEC®plus 翻修股骨头

DE – Deutsch

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Caution! Only valid in conjunction with the instructions for use for ELEC® and ELEC®plus hip joint implants

ELEC®plus revision femoral heads with titanium sleeve

An ELEC® plus revision femoral head with separately supplied titanium sleeve may only be paired with an ELEC®plus acetabular insert or a PE cup approved by the manufacturer of the other components or an approved PE acetabular insert.

The instructions for use from the manufacturer of the other components (endoprosthesis manufacturer) must be observed.

1 Product description and implant raw materials

Please read the corresponding section in the instructions for use for ELEC® and ELEC®plus hip joint implants beforehand.

ELEC®plus revision femoral heads with Ti sleeve are intended to replace a femoral head component by being attached to a stem remaining in situ during revision surgery, which became necessary due to cup loosening or femoral head breakage. ELEC®plus revision femoral heads with Ti sleeve can also be used in primary operations for femoral neck length correction. ELEC®plus revision femoral heads with Ti sleeve should only be used for fixing to stems approved by the endoprosthesis manufacturer.

ELEC®plus revision femoral heads with Ti sleeve are manufactured from high-purity implant ceramic (ISO 6474-2) or a high-strength titanium material (ISO 5832-3).

2 Handling

Please read the corresponding section in the instructions for use for ELEC® and ELEC®plus hip joint implants beforehand.

2.1 General instructions

Caution ELEC®plus revision femoral heads with Ti sleeve are intended for single use only. It is prohibited to reuse a prosthesis component that had previously been implanted in the body. Reuse may lead to breakage due to scratches or impact marks.

2.2 Instructions for use

Only the system instruments recommended by the manufacturer of the other components should be used for the implantation. The enclosed instructions for use should be noted before the instruments are used. ELEC®plus revision femoral heads with Ti sleeve may only be used in the combinations recommended by the manufacturer of the other components. Only ELEC®plus revision femoral heads with Ti sleeve with taper dimensions that match those of the combination implant may be used. The taper size is indicated on the label and, if possible, also on the implant. Only brand new implants should be used.

Assembling ELEC®plus revision femoral heads and Ti sleeve

The ELEC®plus revision femoral heads must be connected with the separately packaged Ti sleeve prior to use. To do this, hold the Ti sleeve at its widest point with two fingers and slide it into the taper of the femoral head with a little rotation until it is firmly positioned. When assembling the components, it is important to ensure that they are completely clean.

2.3 Surgical technique

Attaching the femoral head and Ti sleeve assembly onto the prosthesis stem remaining in situ

In order to guarantee reliable function, it is essential that the utmost care be taken. The following tasks should be completed before attaching the ELEC®plus revision femoral head and Ti sleeve assembly:

1. Detach the femoral head to be replaced with a suitable detaching instrument to avoid

- any damage to the stem taper remaining in situ.
2. Rinse the stem taper remaining in situ thoroughly with water to remove any tissue, bone fragments or cement residue
 3. Inspect the stem taper remaining in situ for damage, if undamaged, continue to 5)
 4. If damaged, assess the damage:
 - a. Permissible damage for the use of an ELEC[®]plus revision femoral head: Fine damage marks on the used stem taper caused when the previous femoral head was detached, minor scratches measuring up to 0.2mm in width and depth
 - b. Unacceptable damage for the use of an ELEC[®]plus revision femoral head: Squashed taper, flattened taper, bevelled taper,
 5. Dry the stem taper with a clean gauze pad and check the position and function of the joint after attaching the test head.
 6. Remove the test head and repeat the rinsing and drying steps

The following tasks should be completed when attaching the ELEC[®]plus revision femoral head and Ti sleeve assembly:

7. Place the ELEC[®]plus revision femoral head and Ti sleeve assembly onto the stem taper with a little axial pressure and simultaneous rotation until it is firmly in place
8. Now place the plastic impactor on the pole of the head and finally clamp the femoral head onto the stem with a light hammer blow in the axial direction. The hammer blow results in the non-twisting fixation of the implant.
9. Check the fixation by trying to remove the implant by hand.

Caution Never hit a ceramic femoral head with a metal hammer.

3 Packaging and sterility

Please read the corresponding section in the instructions for use for ELEC[®] and ELEC[®]plus hip joint implants beforehand.

4 Preoperative planning

Please read the corresponding section in the instructions for use for ELEC[®] and ELEC[®]plus hip joint implants beforehand.

For revision surgery, it must be ensured in particular that the stem taper dimensions of the stem remaining in situ are compatible with the internal taper dimensions of the Ti sleeve of an ELEC[®]plus revision femoral head. The specifications of the endoprosthesis manufacturer must be observed.

5 Indications

- Total hip replacement, primary surgery in combination with prosthesis stems and cups
- Revision (= endoprosthesis removal) using new prosthesis stems and cups, or in situ remaining stems without damage or with permissible damage
- Hemiarthroplasty (partial hip replacement, only applicable to femoral heads)

6 Contraindications

- Revision with retention of stems in situ with unacceptably damaged stem tapers (risk of component breakage)

7 Possible adverse effects and complications

Please read the corresponding section in the instructions for use for ELEC[®] and ELEC[®]plus hip joint implants beforehand.

8 Patient information, documentation

Please read the corresponding section in the instructions for use for ELEC[®] and ELEC[®]plus hip joint implants beforehand.

9 Summary of Safety and Clinical Performance

Please read the corresponding section in the instructions for use for ELEC[®] and ELEC[®]plus hip joint implants beforehand.

10 Explanation of the label symbols:

Please read the corresponding section in the instructions for use for ELEC® and ELEC®plus hip joint implants beforehand.

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