OHST

medical technology

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Grünauer Fenn 3

14712 Rathenow

Deutschland / Germany / Allemagne /

Germania / Alemania / Duitsland /

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IMPLANT

Müller II Cup

Before using the product, the user is under obligation to carefully study the following recommendations and information along with the information specific to the product.

The party introducing this product into circulation accepts no liability for direct or consequential damage or injury resulting from careless use or handling, particularly noncompliance with the following user instructions or improper care or maintenance. These implants may be used only by physicians with appropriate experience and practice in hip arthroplasty. Familiarity with the surgical technique recommended for this system and its diligent application are indispensable in order to achieve the best possible result.

1. Product description and implant materials

The Müller II Cup is a prosthetic hip cup designed for cemented fixation in the acetabulum. It is manufactured using UHMWPE according to ISO 5834-2 and is available in the variants "Standard", "Flat", "Snap" and "Dysplasia". The Standard, Snap and Flat variants are available for spherical head diameters of 22 mm, 28 mm and 32 mm, while the Dysplasia variant is designed to retain femoral heads with a diameter of 28 mm and 32 mm. To ensure radiographic visibility, all implants are equipped with a ring made of implant steel (ISO 5832-1).

Product, packaging contents and materials used are specified on the product label. The implant must be implanted using a suitable surgical technique familiar to the surgeon. In this regard, attention must be given to the explanations concerning the particular surgical technique.

1.1 Implant overview

Name	Material	Reference number
Müller II Cup Standard		
Mueller II Cup (Standard) Ø36, KD22	ISO 5834-2 UHMWPE	1300-22-36
Mueller II Cup (Standard) Ø38, KD22	ISO 5834-2 UHMWPE	1300-22-38
Mueller II Cup (Standard) Ø40, KD22	ISO 5834-2 UHMWPE	1300-22-40
Mueller II Cup (Standard) Ø42, KD22	ISO 5834-2 UHMWPE	1300-22-42
Mueller II Cup (Standard) Ø44, KD22	ISO 5834-2 UHMWPE	1300-22-44
Mueller II Cup (Standard) Ø42, KD28	ISO 5834-2 UHMWPE	1300-28-42
Mueller II Cup (Standard) Ø44, KD28	ISO 5834-2 UHMWPE	1300-28-44
Mueller II Cup (Standard) Ø46, KD28	ISO 5834-2 UHMWPE	1300-28-46
Mueller II Cup (Standard) Ø48, KD28	ISO 5834-2 UHMWPE	1300-28-48
Mueller II Cup (Standard) Ø50, KD28	ISO 5834-2 UHMWPE	1300-28-50
Mueller II Cup (Standard) Ø52, KD28	ISO 5834-2 UHMWPE	1300-28-52
Mueller II Cup (Standard) Ø54, KD28	ISO 5834-2 UHMWPE	1300-28-54
Mueller II Cup (Standard) Ø56, KD28	ISO 5834-2 UHMWPE	1300-28-56

Name	Material	Reference number
Mueller II Cup (Standard) Ø58, KD28	ISO 5834-2 UHMWPE	1300-28-58
Mueller II Cup (Standard) Ø60, KD28	ISO 5834-2 UHMWPE	1300-28-60
Mueller II Cup (Standard) Ø62, KD28	ISO 5834-2 UHMWPE	1300-28-62
Mueller II Cup (Standard) Ø64, KD28	ISO 5834-2 UHMWPE	1300-28-64
Mueller II Cup (Standard) Ø66, KD28	ISO 5834-2 UHMWPE	1300-28-66
Mueller II Cup (Standard) Ø68, KD28	ISO 5834-2 UHMWPE	1300-28-68
Mueller II Cup (Standard) Ø70, KD28	ISO 5834-2 UHMWPE	1300-28-70
Mueller II Cup (Standard) Ø44, KD32	ISO 5834-2 UHMWPE	1300-32-44
Mueller II Cup (Standard) Ø46, KD32	ISO 5834-2 UHMWPE	1300-32-46
Mueller II Cup (Standard) Ø48, KD32	ISO 5834-2 UHMWPE	1300-32-48
Mueller II Cup (Standard) Ø50, KD32	ISO 5834-2 UHMWPE	1300-32-50
Mueller II Cup (Standard) Ø52, KD32	ISO 5834-2 UHMWPE	1300-32-52
Mueller II Cup (Standard) Ø54, KD32	ISO 5834-2 UHMWPE	1300-32-54
Mueller II Cup (Standard) Ø56, KD32	ISO 5834-2 UHMWPE	1300-32-56
Mueller II Cup (Standard) Ø58, KD32	ISO 5834-2 UHMWPE	1300-32-58
Mueller II Cup (Standard) Ø60, KD32	ISO 5834-2 UHMWPE	1300-32-60
Mueller II Cup (Standard) Ø62, KD32	ISO 5834-2 UHMWPE	1300-32-62
Mueller II Cup (Standard) Ø64, KD32	ISO 5834-2 UHMWPE	1300-32-64
Mueller II Cup (Standard) Ø66, KD32	ISO 5834-2 UHMWPE	1300-32-66
Mueller II Cup (Standard) Ø68, KD32	ISO 5834-2 UHMWPE	1300-32-68
Mueller II Cup (Standard) Ø70, KD32	ISO 5834-2 UHMWPE	1300-32-70
Müller II Cup Flat		
Mueller II Cup (Flat) Ø36, KD22	ISO 5834-2 UHMWPE	1301-22-36
Mueller II Cup (Flat) Ø38, KD22	ISO 5834-2 UHMWPE	1301-22-38
Mueller II Cup (Flat) Ø40, KD22	ISO 5834-2 UHMWPE	1301-22-40
Mueller II Cup (Flat) Ø42, KD22	ISO 5834-2 UHMWPE	1301-22-42
Mueller II Cup (Flat) Ø44, KD22	ISO 5834-2 UHMWPE	1301-22-44
Mueller II Cup (Flat) Ø42, KD28	ISO 5834-2 UHMWPE	1301-28-42
Mueller II Cup (Flat) Ø44, KD28	ISO 5834-2 UHMWPE	1301-28-44
Mueller II Cup (Flat) Ø46, KD28	ISO 5834-2 UHMWPE	1301-28-46
Mueller II Cup (Flat) Ø48, KD28	ISO 5834-2 UHMWPE	1301-28-48
Mueller II Cup (Flat) Ø50, KD28	ISO 5834-2 UHMWPE	1301-28-50
Mueller II Cup (Flat) Ø52, KD28	ISO 5834-2 UHMWPE	1301-28-52
Mueller II Cup (Flat) Ø54, KD28	ISO 5834-2 UHMWPE	1301-28-54
Mueller II Cup (Flat) Ø56, KD28	ISO 5834-2 UHMWPE	1301-28-56
Mueller II Cup (Flat) Ø58, KD28	ISO 5834-2 UHMWPE	1301-28-58
Mueller II Cup (Flat) Ø60, KD28	ISO 5834-2 UHMWPE	1301-28-60
Mueller II Cup (Flat) Ø62, KD28	ISO 5834-2 UHMWPE	1301-28-62

Name	Material	Reference number
Mueller II Cup (Flat) Ø64, KD28	ISO 5834-2 UHMWPE	1301-28-64
Mueller II Cup (Flat) Ø66, KD28	ISO 5834-2 UHMWPE	1301-28-66
Mueller II Cup (Flat) Ø68, KD28	ISO 5834-2 UHMWPE	1301-28-68
Mueller II Cup (Flat) Ø70, KD28	ISO 5834-2 UHMWPE	1301-28-70
Mueller II Cup (Flat) Ø44, KD32	ISO 5834-2 UHMWPE	1301-32-44
Mueller II Cup (Flat) Ø46, KD32	ISO 5834-2 UHMWPE	1301-32-46
Mueller II Cup (Flat) Ø48, KD32	ISO 5834-2 UHMWPE	1301-32-48
Mueller II Cup (Flat) Ø50, KD32	ISO 5834-2 UHMWPE	1301-32-50
Mueller II Cup (Flat) Ø52, KD32	ISO 5834-2 UHMWPE	1301-32-52
Mueller II Cup (Flat) Ø54, KD32	ISO 5834-2 UHMWPE	1301-32-54
Mueller II Cup (Flat) Ø56, KD32	ISO 5834-2 UHMWPE	1301-32-56
Mueller II Cup (Flat) Ø58, KD32	ISO 5834-2 UHMWPE	1301-32-58
Mueller II Cup (Flat) Ø60, KD32	ISO 5834-2 UHMWPE	1301-32-60
Mueller II Cup (Flat) Ø62, KD32	ISO 5834-2 UHMWPE	1301-32-62
Mueller II Cup (Flat) Ø64, KD32	ISO 5834-2 UHMWPE	1301-32-64
Mueller II Cup (Flat) Ø66, KD32	ISO 5834-2 UHMWPE	1301-32-66
Mueller II Cup (Flat) Ø68, KD32	ISO 5834-2 UHMWPE	1301-32-68
Mueller II Cup (Flat) Ø70, KD32	ISO 5834-2 UHMWPE	1301-32-70
Müller II Cup Snap		
Mueller II Cup (Snap) Ø36, KD22	ISO 5834-2 UHMWPE	1310-22-36
Mueller II Cup (Snap) Ø38, KD22	ISO 5834-2 UHMWPE	1310-22-38
Mueller II Cup (Snap) Ø40, KD22	ISO 5834-2 UHMWPE	1310-22-40
Mueller II Cup (Snap) Ø42, KD22	ISO 5834-2 UHMWPE	1310-22-42
Mueller II Cup (Snap) Ø44, KD22	ISO 5834-2 UHMWPE	1310-22-44
Mueller II Cup (Snap) Ø44, KD28	ISO 5834-2 UHMWPE	1310-28-42
Mueller II Cup (Snap) Ø42, KD28	ISO 5834-2 UHMWPE	1310-28-44
Mueller II Cup (Snap) Ø46, KD28	ISO 5834-2 UHMWPE	1310-28-46
Mueller II Cup (Snap) Ø48, KD28	ISO 5834-2 UHMWPE	1310-28-48
Mueller II Cup (Snap) Ø50, KD28	ISO 5834-2 UHMWPE	1310-28-50
Mueller II Cup (Snap) Ø52, KD28	ISO 5834-2 UHMWPE	1310-28-52
Mueller II Cup (Snap) Ø54, KD28	ISO 5834-2 UHMWPE	1310-28-54
Mueller II Cup (Snap) Ø56, KD28	ISO 5834-2 UHMWPE	1310-28-56
Mueller II Cup (Snap) Ø58, KD28	ISO 5834-2 UHMWPE	1310-28-58
Mueller II Cup (Snap) Ø60, KD28	ISO 5834-2 UHMWPE	1310-28-60
Mueller II Cup (Snap) Ø62, KD28	ISO 5834-2 UHMWPE	1310-28-62
Mueller II Cup (Snap) Ø64, KD28	ISO 5834-2 UHMWPE	1310-28-64
Mueller II Cup (Snap) Ø66, KD28	ISO 5834-2 UHMWPE	1310-28-66
Mueller II Cup (Snap) Ø68, KD28	ISO 5834-2 UHMWPE	1310-28-68

Name	Material	Reference number
Mueller II Cup (Snap) Ø70, KD28	ISO 5834-2 UHMWPE	1310-28-70
Mueller II Cup (Snap) Ø44, KD32	ISO 5834-2 UHMWPE	1310-32-44
Mueller II Cup (Snap) Ø46, KD32	ISO 5834-2 UHMWPE	1310-32-46
Mueller II Cup (Snap) Ø48, KD32	ISO 5834-2 UHMWPE	1310-32-48
Mueller II Cup (Snap) Ø50, KD32	ISO 5834-2 UHMWPE	1310-32-50
Mueller II Cup (Snap) Ø52, KD32	ISO 5834-2 UHMWPE	1310-32-52
Mueller II Cup (Snap) Ø54, KD32	ISO 5834-2 UHMWPE	1310-32-54
Mueller II Cup (Snap) Ø56, KD32	ISO 5834-2 UHMWPE	1310-32-56
Mueller II Cup (Snap) Ø58, KD32	ISO 5834-2 UHMWPE	1310-32-58
Mueller II Cup (Snap) Ø60, KD32	ISO 5834-2 UHMWPE	1310-32-60
Mueller II Cup (Snap) Ø62, KD32	ISO 5834-2 UHMWPE	1310-32-62
Mueller II Cup (Snap) Ø64, KD32	ISO 5834-2 UHMWPE	1310-32-64
Mueller II Cup (Snap) Ø66, KD32	ISO 5834-2 UHMWPE	1310-32-66
Mueller II Cup (Snap) Ø68, KD32	ISO 5834-2 UHMWPE	1310-32-68
Mueller II Cup (Snap) Ø70, KD32	ISO 5834-2 UHMWPE	1310-32-70
Müller II Cup Dysplasia		
Mueller II Cup (Dysplasia) Ø42, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-42
Mueller II Cup (Dysplasia) Ø44, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-44
Mueller II Cup (Dysplasia) Ø46, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-46
Mueller II Cup (Dysplasia) Ø48, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-48
Mueller II Cup (Dysplasia) Ø50, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-50
Mueller II Cup (Dysplasia) Ø52, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-52
Mueller II Cup (Dysplasia) Ø54, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-54
Mueller II Cup (Dysplasia) Ø56, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-56
Mueller II Cup (Dysplasia) Ø58, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-58
Mueller II Cup (Dysplasia) Ø60, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-60
Mueller II Cup (Dysplasia) Ø62, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-62
Mueller II Cup (Dysplasia) Ø64, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-64
Mueller II Cup (Dysplasia) Ø66, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-66
Mueller II Cup (Dysplasia) Ø68, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-68
Mueller II Cup (Dysplasia) Ø70, KD28 - 10°	ISO 5834-2 UHMWPE	1303-28-70
Mueller II Cup (Dysplasia) Ø44, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-44
Mueller II Cup (Dysplasia) Ø46, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-46
Mueller II Cup (Dysplasia) Ø48, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-48
Mueller II Cup (Dysplasia) Ø50, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-50
Mueller II Cup (Dysplasia) Ø52, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-52
Mueller II Cup (Dysplasia) Ø54, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-54
Mueller II Cup (Dysplasia) Ø56, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-56

Name	Material	Reference number
Mueller II Cup (Dysplasia) Ø58, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-58
Mueller II Cup (Dysplasia) Ø60, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-60
Mueller II Cup (Dysplasia) Ø62, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-62
Mueller II Cup (Dysplasia) Ø64, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-64
Mueller II Cup (Dysplasia) Ø66, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-66
Mueller II Cup (Dysplasia) Ø68, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-68
Mueller II Cup (Dysplasia) Ø70, KD32 - 10°	ISO 5834-2 UHMWPE	1303-32-70

1.2 Instrument overview

The instruments of the OHST Medizintechnik AG listed below must be used exclusively for implantation:

Name	Reference number
Acetabulum Reamer Instrument Set	367-147
Equatorial Size Testers Instrument Set	367-148
Mueller II Cup Insertion Instrument Set	367-400
Optional:	
Mueller II Cup Insertion Instrument Set (small)	367-402
Cup Positioning Guide Ø8 to Ø24	506-010
Acetabular Reamer Ø40 mm	506-540

1.3 Accessories

Name	Reference number
Surgical Technique for Müller II Cup	50000402
X-ray Templates Mueller II Cup Dysplasia	367-2001
X-ray Templates Mueller II Cup Standard KD 22 mm	367-2002
X-ray Templates Mueller II Cup Standard	367-2003
X-ray Templates Mueller II Cup Flat	367-2004
X-ray Templates Mueller II Cup Snap	367-2005
X-ray Templates Mueller II Cup Snap KD22 mm	367-2006
X-ray Templates Mueller II Cup Flat KD 22 mm	367-2010
Implant passport	50000572

2. Handling

2.1 General information

This implant is part of a system and must only be used with the appropriate original system components. Only the instruments of the system listed above must be used for implantation. Before using the instruments the attached instructions for use (50000354) must be considered.

Caution: Implants must always be kept in their complete, unopened protective packaging. The packaging containing the implant must not be exposed to direct sunlight. Before inserting the implant, the packaging must be examined for damage, as this could affect sterility.

When unpacking the implant, its conformity with the designation on the packaging (art. no. / serial no. / size) must be checked. Compliance is required with the appropriate hygiene regulations during removal of the implant from the packaging. Care must be taken to protect all implant surfaces against damage, since this could be decisive for possible failure. The prosthesis must not therefore come into contact with objects which could damage its surface. Before use, every implant must be visually inspected for damage. Machining an implant can not only reduce its service life, but can also lead to immediate or subsequent failure of the prosthesis under stress. The implant must therefore neither be mechanically nor otherwise processed. Implants from damaged packaging, unsterile, contaminated, damaged or carelessly handled implants or implants subjected to unauthorized machining must not be used.

Caution: Implants are intended for single use only! Individual loads on functional surfaces of an implant used for one patient modify the functional surfaces in a way that excludes any reuse. Detection of load-caused markings by visual methods only is not secured. Therefore, damage after explantation must be assumed which excludes any reuse.

2.2 Authorised component combinations

We only guarantee compatibility of our products in combination with our own CE-marked products and with the products we have approved for combined use, for which the competent authority has issued its approval. In this regard, please note the instructions for use of the endoprosthesis manufacturers and the combination matrix approved by OHST.

Due to reasons relating to the product safety and product liability, it is prohibited to use implants manufactured by OHST Medizintechnik AG in combination with components of other manufacturer that have not been approved by OHST.

2.3 Information for use

The application of the implant is done with cement. For the proper accomplishment of the cementing step the manufacturer data about the application of the bone cement must be considered. To lower the risk of severe cardiovascular complications (caused by BCIS = bone cement implantation syndrome), we recommend using vacuum-mixed bone cement.

Caution: The ball diameter of the femoral head prosthesis must agree perfectly with the spherical nominal diameter of the part of the acetabular cup implant which articulates as the sliding partner.

Caution: Please be hereby explicitly advised that, in case of an intraoperative change or revision of the femoral head, only femoral heads without a ceramic cone are to be used. This is valid irrespective of the materials used in the previous cone pairing.

Caution: If a ceramic component is damaged or fractured, complete revision of the prosthetic components at the earliest possible date is recommended. In this case, the use of metal femoral heads is contraindicated in revision surgery, as this may lead to serious and partly life-threatening complications. In the rare event of a fracture of the ceramic component, thorough debridement with removal of all visible ceramic particles as well as careful wound irrigation is absolutely essential during surgery.

Before inserting the implant, the implant bed must be irrigated. During implantation, ensure that all loose particles (e. g. bone splinters, friction particles from the instruments) are removed from the prepared implant bed.

Caution: When using high-frequency surgical instruments (e.g. cautery knife), it must be ensured that they do not come into contact with the implants or instruments. This can cause such severe damage to the implants or instruments that failure (e.g. fracture) may ensue. If an implant has been damaged, it must not remain in the patient but needs to be replaced by a new, intact implant. Damaged instruments may only continue to be used if they can still perform their intended function without compromise.

2.4 Surgical technique

Any approach to the hip joint deemed suitable by the surgeon can be used. The surgeon should have a good view of the anatomical structures so that correct working with the instruments is not impeded. After opening the joint capsule and dislocating the femoral head from the acetabulum, resect the femoral head as determined in the preoperative planning of the hip stem.

Starting with the smallest, ream the acetabulum using spherical acetabulum reamers in ascending order; in doing so, make sure to always let the reamers move freely and never apply force. Always allow the reamer to move freely and do not distort the reamer sizes by applying pressure laterally! The base of the cup should be completely free from cartilage and the subchondral bone should bleed evenly. The bony mass collected from the last reaming process can be used to fill in any potential gaps between the implant and the acetabulum.

Caution: The diameter (nominal size) of the last acetabulum reamer used has to be at least 4 mm larger than the smallest Müller II Cup available to ensure a cement coating of 2 mm all around.

Next, insert the respective equatorial size tester based on the size of the reamer and use it to check the alignment (45° inclination and 10° anteversion) of the hip cup. The equatorial size tester should remain in place through suction and be surrounded by a sufficient quantity of bone. After these surgical steps are complete, select the appropriate implant size.

Follow the cement manufacturer's instructions for use to prepare the bone cement using a state-of-the-art cementing technique and insert it into the prepared implant bed.

Caution: Make sure that a cement coating of around 2 mm is achieved. Use utmost care when removing excess cement to prevent damaging the cup or to prevent cement particles from penetrating the inner articulation surface.

The size of the implant to be used has to be 4 mm smaller than the size of the equatorial size tester selected to ensure an adequate cement coating. Implant the Müller II Cup with the corresponding insertion head that is screwed on to the curved handle. To this purpose, press the implant into the bone cement and align it in accordance with the pre-operative planning process. Then apply light pressure and wait until the bone cement has set. This completes the insertion of the cup. In order to protect the cup against damage during the remainder of the operation it is advisable to cover it up with a sterile compress. After all of the other components have been implanted, close the wound using the routine procedure.

When using the snap variant of the Müller II cup, greater resistance has to be overcome to reduce the femoral head into the cup than with the other Müller II cup variants. When using the snap variant of the Müller II cup, the femoral head can only be reduced after flexion/anteversion of the thigh until the chamfer in the neck area of the femoral head protrudes into the internal contour of the cup, allowing air to escape from the joint body. Otherwise an air pocket would form, thus preventing the femoral head from snapping into place. After reducing the femoral head, verify the correct fit of the femoral head inside the cup. The connection between the two implants is correct when the femoral head noticeably snaps into the cup.

3. Packaging and sterility

Depending on the sterilisation method used, implants are packaged in a triple transparent pouch made of plastic laminated film (sterilisation by irradiation at least 25 kGy) or in a double transparent pouch made of Tyvek® (ethylene oxide sterilisation) with a carton. The instruments are supplied unsterile in protective packaging. They must be cleaned and sterilised prior to use in accordance with the respective instructions for use (50000354). The stated expiry date presumes that the packaging is intact and unopened and that the product is stored under suitable conditions.

Caution: The implants may not be resterilized! The reconditioning of components that have not been implanted but the packaging of which has been opened is permitted only at the manufacturer, because the components must pass through individual validated processes once again.

The outer pouch of the triple transparent pouch packaging is to be removed by the non-sterile personnel together with the carton. For the double transparent pouch packaging, only the carton is to be removed by the non-sterile personnel. The second pouch must be opened such that the sterility of the inner pouch is not compromised. The inner pouch is removed and opened by the sterile personnel. The implant must then be presented to the surgeon, who can then directly remove the sterile implant.

4. Preoperative planning and postoperative care

Preoperative planning by reference to X-rays, CT data and similar is indispensable and provides important information about suitable implants, placing and possible component combinations and enables the size of the implant to be used to be preselected. Surgery may only be performed once it has been established that the patient is able to tolerate the implant material. Use the X-ray templates for planning the operation. These are available for all sizes in a magnification of 1,15:1. In addition, X-ray templates with a 1:1 ratio are available in digital form. Trial prostheses for checking the correct seating (where applicable) and additional implants should be available should another size be required or the intended implant cannot be used. Recognized procedures must be used for postoperative care.

5. Indications

- Advanced degeneration of the hip joint due to degenerative, post-traumatic or rheumatoid arthritis
- Avascular necrosis of the femoral head
- Displaced intracapsular hip fracture
- Sequelae of earlier surgical procedures (osteosynthesis, articular reconstruction, arthrodesis, etc.)

6. Contraindications

- Acute local (in relation to the joint undergoing surgery) or active systemic infection
- Insufficient bone substance or quality that does not allow for stable fixation of the implant
- Hypersensitivity to the materials used
- Severe osteoporosis
- Severe malformations, congenital hip dislocations
- Local bone tumours that preclude stable fixation of the implant
- Growth in children and adolescents

When using Snap or Dysplasia cups or cup inserts, the movement range in flexion and extension is reduced by approximately 55° (Snap) and 27° (Dysplasia) and in abduction and adduction by approximately 22° (Snap) and 10° (Dysplasia) as compared to the standard application.

7. Risk factors and conditions that may affect the success of the surgery

Caution: Clinical experience has shown that the presence of one or more of the following concomitant circumstances (risk factors) may lead to shorter service lives, more frequent complications or an altogether poorer outcome of hip arthroplasty. This list is by no means exhaustive.

General risk factors and conditions:

- Overweight
- Alcohol or substance abuse
- Patient groups with mental disorders or addictions
- Pregnancy
- High-dose ingestion of cortisone or cytostatics
- Previous or threatening infectious diseases with possible joint involvement
- Deep vein thrombosis and/or history of pulmonary embolism
- Systemic diseases and metabolic disorders
- All general surgical risks

Risk factors and conditions specific to hip arthroplasty:

- Occurrence of fissures, in rare cases fractures
- Circulatory disorders of the affected limb
- Neurological disorders of the affected limb
- Muscle malfunction in the affected limb
- Muscle spasms or other spastic conditions
- Anticipated extreme loading e.g. due to work and sport
- Epilepsy or other reasons for repeated trauma with an increased risk of fracture
- Joint deformities that make fixation of the implant difficult
- Weakening of the bearing structures by tumour (e.g. bone cyst, non-ossifying fibroma)
- Loss of the ligamenture

8. Possible negative effects

The negative effects listed below are among the most typical and commonly occurring consequences of a surgical procedure:

- Infection
- Venous thrombosis and pulmonary embolism
- Cardiovascular disorders
- Haematomas
- Paresthesia
- Numbness
- Swelling
- Nerve damage
- Oedem/fluids

The negative effects listed below are among the most typical and most frequently occurring consequences of total hip arthroplasty:

- Change in position and loosening of the prosthesis
- Dislocation of the prosthesis
- Implant breakage
- Stiffness
- Reduced quality of life (pain, sleep disorders, ROM limitations; in particular also when lying down)
- Osteolysis
- Heterotopic ossification
- Cardiac arrhythmia, increased pulmonary vascular resistance, cardiac arrest (caused by BCIS = bone cement implantation syndrome)

Caution: Revision surgery may be required if specific undesirable effects occur.

9. Patient information, documentation

The serial numbers of the implants used must be recorded in the patient's records. Appropriate labels are included with the packaging of the sterile implants.

The patient must be informed of the advantages and risks of the procedure. If the implant is regarded as the best solution for the patient, although the contraindications described above partially apply to the patient, it is particularly important to point out to the patient the effects of these circumstances on the success. It must be explained to patients receiving a hip replacement that the life of the implant will depend on their weight and degree of activity. The patient has to be informed about activities with which he can reduce the effects of these aggravating circumstances.

All the information given to the patient must be documented in writing by the surgeon. After surgery, the patient must be given an implant pass containing all necessary information concerning the implant. Adhesive labels are enclosed for documenting the implant used. Adverse effects that are harmful to patients can arise during MRI investigations. Artefacts, heating of implant, induction of electrical currents and implant loosening are among the possible effects. The equipment manufacturer's instructions should be carefully studied before use. In case of doubt, comparable implants should be checked for their specific MRI suitability as part of an individual risk assessment. Patients should be informed of the risks.

10. Key to label symbols

An explanation of the symbols used by OHST Medizintechnik AG can be found in the annex (p. 217).

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