CEMSTOP®

Obturateur diaphysaire biorésorbable
Bio-absorbable cement restrictor
Biore sorbierbare Markraumstopper
Obturador diafisario bioresorbible
Otturatore Diafisario biorassorbibile
Restrictor de Cimento bio-absorbível
Biyorezorbabil diafiz kanalı obturatori
Enxerto para Ossos e Dispositivos associados

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Before using TEKNIMED products, the operator must thoroughly familiarise himself with the safety regulations figuring in the directions for use, as well as the information relating to each product (description, operating technique, leaflets, etc.). The corresponding information is available from TEKNIMED. The surgeon should also be informed of the residual risk associated with the product he considers using.

GENERAL INFORMATION

Implantation of TEKNIMED products should only be performed by qualified operators having a sound knowledge and full mastery of operative techniques specific to TEKNIMED products. Operative techniques may be acquired from the TEKNIMED distributors. The surgeon is responsible for any complications or harmful consequences which might result from an erroneous indication or operative technique, improper use of the equipment and failure to comply with the safety regulations provided in the directions for use. The cement stopper is delivered sterile and ready for use.

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RECOMMANDATIONS FOR DISPOSAL

The device and its components must be discarded in compliance with local regulations on waste disposal.

SÉCURITÉ

The diameter of the CEMSTOP® can be pre-determined by the surgeon by setting up the correspondence between the supplied templates and diameters of the medullary canal.

COMPOSITION

Gelatine (porcine-based)

Glycerol

Water

Methylparahydroxybenzoate

These components comply with European Pharmacopoeia standards.

INSTRUCTIONS FOR USE

The diameter of the CEMSTOP® can be pre-determined by the surgeon by setting up the correspondence between the supplied templates and diameters of the medullary canal.

STERILISATION

CEMSTOP® is sterilized by gamma radiation at a minimum dose of 25 kGy.

Prior to use, carefully check the protective packaging in order to ensure that it has not been damaged in a way that could affect its sterility.

When taking the product out of its packaging, be aware of the asepsis rules.

CEMSTOP® is delivered sterile and ready for use.

The surgeon is responsible for any complications or harmful consequences which might result from an erroneous indication or operative technique, improper use of the equipment and failure to comply with the safety regulations provided in the directions for use.

CEMSTOP® is completely reabsorbed within several days at 37°C

and resists the temporary increase in temperature generated by the polymerization of acrylic cement.

PATIENT INFORMATION

The patient should be informed by the doctor of the potential consequences of the factors mentioned in the following paragraphs: contra-indications and side effects, that is, those liable to hinder the success of the operation, as well as possible complications which may arise. The patient should also be informed of the measures to be taken to diminish the possible consequences of these factors.

INDICATIONS

CEMSTOP® cement restrictor is a diaphyseal plug designed to occlude the medullary cavity before the introduction of acrylic cement during total hip arthroplasty.

CEMSTOP® prevents the cement from flowing down the diaphysis, and therefore facilitates cement pressurization. Completely biocompatible, CEMSTOP® is reabsorbed in a plasto

to several days at 37°C and resists the temporary increase in temperature generated by the polymerization of acrylic cement.

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