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Option Kugelkopf

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Lit. No. 848001 Rev. 0 As from 03/2018

**1. Product Name**

The Option ball head consists of a shell (material: titanium alloy TiAl6V4 in accordance with ISO 5832-3 and ASTM F136) with inner cone 12/14 and a BIOLOX® delta ceramic ball head with a conic drill-hole to match the shell.

**2. Indications and Contraindications**

The Option ball head is used in order to implant a ceramic ball head again in cup revisions or when replacing a ball head on a prosthetic stem from the company ImplanTec GmbH that is left in situ.

The Option ball head can be used in the extremely rare event of a ceramic component breaking.

The Option ball head can also be used in primary surgery for neck length corrections.

The indications and contraindications must comply with the general rule and the current state of the art for hip replacement surgery. The recommendations and notes on indications and contraindications in the instructions for use for the ANA.NOVA® hip system, Lit. No. 841001, must be observed.

**3. Properties of BIOLOX® delta Ceramic**

BIOLOX® delta is an ultra-pure mixed aluminium oxide ceramic in accordance with ISO 6474-2. It is bio-inert, bio-compatible, bio-stable, mechanically stable, corrosion-proof, does not cause any allergic reactions, and displays excellent resistance to fatigue, a high degree of impact resistance and tensile strength, excellent fracture toughness and an extremely high degree of hardness. A further - clinically proven - major advantage of BIOLOX® delta ball heads is the fact that by comparison with metal ball heads there is less polyethylene abrasion when combined with polyethylene cups or cup inserts.

**4. Important information for the Surgeon**

The Option ball head was sterilised with gamma radiation of at least 25 kGy, as indicated by the labelling on the packaging and the discoloured sterilisation indicator.

An Option ball head that is delivered sterile must be kept unopened in the original packaging until use. Before use, the sterile expiry date on the label must be observed and the protective packaging checked for damage.

The Option ball head must not be used if the sterile expiry date has expired or the protective packaging is damaged.

The colour of an Option ball head sterilized by gamma radiation may be altered. This has no influence on the strength or any other properties of BIOLOX® delta ceramic.

The Option ball head must not be re-sterilised. In exceptional cases, re-sterilisation of the Option ball head may be performed only by ImplanTec GmbH, because only they can comply with the process steps and parameters that are necessary for achieving sterility successfully, and for the prevention of possible impairments of the ceramic.

**Reusability and Prevention of Damage**

Only unused Option ball heads and shells taken from the original packaging may be used. A used Option ball head must not be reused. This means, for instance, that an Option ball head component that has been placed on a stem and then removed again must not be placed on a stem again. An Option ball head that is damaged in any way must not be used either and must be disposed of. This also applies, for instance, to an Option ball head component that has fallen on the floor.

In the case of already used Option ball heads, there is a risk of damage that might not be visible to the naked eye. Damage of any kind can have a negative influence on the functionality and/or stability of the ceramic.

Thus, safe use cannot be guaranteed and therefore only originally packaged, unused and undamaged Option ball head components may be implanted.

**Admissible Combination of Prosthetic Components**

The implants and instruments manufactured and distributed by ImplanTec GmbH for hip replacement can be combined with each other.

Option ball heads may only be combined with the ANA.NOVA® prosthetic components of the company ImplanTec GmbH. Thereby the instructions for use, Lit. No. 841001, must be observed.

On request, ImplanTec GmbH can check and assess the compatibility of a combination with appropriate products from third-party suppliers.

**Assembly of the Option Ball Head**

The Option ball head is delivered as separate components (ceramic ball head, shell). The Option ball head and the shell must be implanted together. Before finally placing the ceramic ball head, the surgeon must assemble the Option ball head system in the packaging tray as shown in Figure 1.

The ceramic ball head is mounted on the shell, which remains in the defined position, by placing it on the shell and pressing it down until resistance is encountered (Fig. 1.1 to 1.4). Attention must be paid not to tilt or twist the ceramic ball head when placing it on the shell.

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
Fig. 1.1 to 1.4 Assembly of the Option ball head

**Preoperative Planning**

The operation is planned based on the available information. Essentially, the stem remaining in situ and the stem cone must be determined during the preoperative planning phase. The inner cone of the shell must correspond with the stem cone. Technical information can be found in the instructions for use of the ANA.NOVA® hip system, Lit. No. 841001.

**Intra-operative Procedure**

**Removing the head and checking the used stem cone**

In revision surgery, the remaining ball head (and shell, if applicable) must be removed using a suitable extraction instrument to avoid any unnecessary damage to the stem cone. After extraction, the remaining stem cone must be checked. If the cone is undamaged, the Option ball head can be used with the cone. If the cone is damaged, the surgeon must check the extent of the damage and satisfy himself that the damage is tolerable.

**Checking the stem cone and decision-making criteria:**

- Permissible condition: used stem cones with minute defects as a result of head-stem disassembly.
- Inadmissible cone deformations: see Fig. II.1 to II.3. Under these circumstances the Option ball head system may not be used with the stem cones.


Fig. II.1 Conical cone, Fig. II.2 Cone with wide flattening, Fig. II.3 Squashed cone

Trials heads must be used to determine the neck length and to check the soft tissue balance and scope of movement.

**Final placement on the stem**

Before final assembly, the stem cone, shell and ceramic ball head must be dry and free from any extraneous matter (e.g. tissue, bone or cement particles). The Option ball head is placed on the stem cone manually with a rotating movement and pressure until it locks.

Generally, it should be possible to place the Option ball head on the stem cone easily. The system must not be used in cases in which the Option ball head can only be placed on the stem cone with pressure.

The surgeon must position the ball head impactor on the pole of the Option ball head and then fix it on the stem cone definitively with a moderate tap on the ball head impactor in axial direction. Multiple moderate taps may also be applied in order to fix the Option ball head and the shell on the stem cone. The surface structure of the metal cone and the stem cone undergoes plastic deformation as a result of tapping with the ball head impactor, resulting in an optimal distribution of pressure and twist-proof fixation.

Successful fixation of the head must be checked by attempting to remove it by hand.

**Caution:** Never tap on the Option ball head with a metal hammer. Only the ball head impactor supplied by ImplanTec GmbH may be used.

**5. Labelling of the Available Option Ball Head Types**

The Option ball head consists of a BIOLOX® delta ball head with diameters from 28 to 48 mm and a metal shell for one of the following dimensions:

S = - 3,0 mm  
M = 0 mm  
L = 4,0 mm  
XL = 7,0 mm

**6. Interactions with Medicinal Products**

Interactions with medicinal products are not known.

**7. Warnings**

The Option ball head may only be used by qualified surgeons with a profound knowledge and experience in the area of hip replacement. The surgeon must observe the instructions for use of the ANA.NOVA® hip system, Lit. No. 841001, before use of the Option ball head and during the operation.

The shell and the ceramic ball head are components of the Option ball head. Only unused Option ball heads may be used. Reuse and re-sterilisation of these components is not permitted.

The surgeon must check the cone of the stem remaining in situ for compatibility with the cone angle of the shell. If there is any doubt, the manufacturer of the stem must be contacted.

In the event of breakage of a ceramic component, a combination of metal (ball head) with polyethylene (insert) as well as metal with metal may not be used for revision surgery.

**8. Patient Behaviour**

The patient must be informed that it is not possible to make the same demands on an artificial joint as on a natural one. Any form of competitive or extreme sports as well as sports involving jerky and sudden movements are contraindicated. The patient must be informed about possible postoperative complications; this information must comply with the current state of medical findings. Patients with unrealistic expectations with regard to the artificial joint, or patients with excess body weight or a weak skeleton who are physically very active are at greater risk.

**9. Pictograms and Figures**

The notes on the packaging and pictograms must be observed.

	Not for reuse.
	Do not re-sterilise
	Batch designation
	Order number
	Sterilisation by radiation
	Do not use if packaging is damaged
	Caution
	Expiry date
	Manufacturer
	Temperature range for storage and operating environment +15°C to +35°C
	Store in a dry place.
	Protect from sunlight
	Observe instructions for use

**10. Information about the Risks of Interfering Interaction**

Undesired effects that harm the patient can occur when performing MRI / CT scans. Possible effects include artefacts, warming of the implant, induction of electrical currents, loosening of the implant. The manufacturer's instructions for use must be studied before the device is used. In case of doubt, reference implants must be tested for suitability in the respective MRI / CT instrument within the scope of an individual risk assessment. The patient must be informed about the risks.

**11. Cleaning, Disinfection and Sterilisation / Re-sterilisation of Instruments**

For information about the treatment of the instruments, see the Cleaning and Treatment Instructions Lit. No. 865001. The manufacturer cannot be held liable for instrument cleaning and sterilisation processes carried out at the buyer's facilities.

**Important:**

For further information about the use of this product please contact your local sales office or visit the homepage [www.implan-tec.at](http://www.implan-tec.at)

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