1. Allgemeine Information

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1. General Information

A2® hip stem system, ELEC®plus ceramic head and ELEC®plus revision ceramic head with titanium sleeve

Before using these products, the user is obliged to carefully study the following recommendations and notes as well as the product-specific instructions.

The person placing these products on the market does not accept any liability for direct or consequential damage resulting from improper use or handling, in particular from noncompliance to the following instructions for use or from improper care or maintenance.
These implants may only be used by surgeons with appropriate experience and practice in hip endoprosthetics.

Product description and package contents

- The product, the package contents and the materials used are specified on the product label. This product must be used in accordance with the relevant product-specific surgical technique.
- The batch and SN number(s) of the implants used must be documented in the patient's records. The corresponding labels for this purpose are enclosed in the packaging of the sterile implants.
- Cementless models are available in varying sizes with a G and a B variant and cemented models are available in varying sizes with a B variant.
- The implant is delivered without bone cement.

Instructions on reuse

- It is not permitted to reuse an implant which has previously been implanted in the body.
- Damage such as scratches, nicks or impact marks can lead to breakage.
- The implants are supplied sterile (sterilised by gamma radiation). They must not be resterilised.

Warnings / precautions

- Warnings on the packaging must be strictly observed.
- This product may only be used in its original condition. Modifications of any nature as well as mechanical manipulations are not permissible.
- The stem models intended for cemented anchorage are labeled "cemeted" and may only be implanted in conjunction with suitable bone cement.
- The stem models intended for cementless anchoring are labeled "non-cemented" and can be recognized by the rough surface structure in the anchoring area.
- The choice whether cemented or cementless is up to the surgeon, whereby the indications, contraindications, surgical techniques and the current guidelines for implanting primary hip stems must be observed.

2. Storage and Handling

Storage of the sterile implants

Implants must always be stored in their unopened original packaging. The implants should be stored in a dry and clean place, protected from direct sunlight at a storage temperature between 10° C and 30° C.
The materials were selected such that short-term temperature fluctuations, which can occur during transportation, do not have a negative impact on the performance of the implants and their packaging.

Handling of the sterile implants

 Prior to insertion of the sterile implant, the packaging should be inspected for damage. If the primary packaging is damaged, the sterility is no longer guaranteed and the product must not be used.

The implants are vacuum-packed in a triple peel-bag packaging (sterilization barrier system). If it is recognizable that air has penetrated into the sterilization barrier system (air-drawn air), the sterility is also no longer guaranteed - the product is not allowed to use.

When the implant is unwrapped, a check is to be carried out to ensure that it corresponds to the designation on the packaging (article no. and size).

The appropriate aseptic regulations for surgical personnel must be observed during removal of the implant from the packaging: The outermost and, if applicable, the middle PE bag are opened by the non-aseptic surgical assistant; the inner bag is opened by the aseptic surgical assistant. The implant must not come into contact with any objects which could damage its surface. Prior to insertion, every implant must be visually inspected for damage.

3. Instructions for Use - Special Application

Instructions Information of the patient

Prior to the surgical procedure, each patient must be informed in detail about the procedure and about possible risks.

A patient receiving a hip replacement should be informed that the life of the implant will depend on the patient's weight and how active the patient is in everyday life. Strenuous physical activity or trauma can lead to loosening, disproportionate wear and/or breakage of the implant. The patient must be informed about the capabilities of the implant and its effects on his/her lifestyle. There is a possibility that the implant will not last for the rest of the patient's life or only for a certain period of time. Implants are not as durable and stable as natural, healthy tissue or natural, healthy bones. It is therefore possible that all components might have to be replaced after a period of time.

Preoperative planning

The preoperative planning provides important information about the suitable type of component, placement and possible component combinations. Additional implants should be kept handy in case other sizes are required or the intended implant cannot be used.

Surgical technique

The implantation must be performed according to the corresponding surgical technique in the valid version. This current version can be downloaded from the ARTIQO GmbH homepage (www.artiqo.de/downloads). Familiarity with the surgical technique recommended for this system and its careful application are essential in order to achieve the best possible results.

Restriction regarding combinations

 Only implants that are supplied by ARTIQO GmbH and which have been approved for combined use can be combined with one an- other for hip replacement procedures. Unsuitable material or product pairings can lead to breakage, premature wear, loosening of the prosthesis, contact corrosion and other consequences.

The surgeon must always ensure that the individual implant components are compatible with one another and observe the general restrictions concerning the combination of materials.

Stainless steel ball heads (ISO 5832-1, ISO 5832-9) may only be combined with titanium stems (ISO 5832-3) to a limited extent in patients with low level of activity and/or advanced age.

Stems made of forged stainless steel (ISO 5832-9) must not be combined with ball heads of cobalt-chromium-molybdenum alloy (ISO 5832-4).

The components of the acetabular joint part (acetabulum / inlay) and of the femoral joint part (stem / femoral head) should always be procured from one supplier.

All hip stems of the company ARTIQO GmbH are manufactured with a stem cone 12/14. The combination with our offered ceramic heads of Biolox®-Delta and Biolox®-Delta Option as well as ELEC®, ELEC®PLUS and ELEC®PLUS Revision (ISO 6474) S/M/L/XL was specially tested.

If separate handling instructions / instructions for use are supplied with the femoral heads, these must also be considered.

On request, the combination with corresponding foreign products can be reviewed by ARTIQO GmbH and the compatibility be assessed.

When using bone cement, the instructions for use of the cement manufacturer must be observed.

 CAUTION: Femoral heads of size XXL or larger are not suitable for combined use with the A2® hip stem due to the increased leverage.

Tested and approved implant combinations:

 **The A2® stem can be combined with the following ball heads, whereby only implants which have been delivered by ARTIQO GmbH and approved for the respective system combination may be used.**

Used materials:

A2 stem cementless (cone): Ti-6Al-4V (ISO 5832-3)

A2 stem cemented (cone): Steel (ISO 5832-9)

The ceramic heads consist of Al₂O₃ or mixed ceramic Al₂O₃ / ZrO₃ (ELEC according to ISO 6474-1, Biolox®Delta and Biolox® option or ELEC®Plus and ELEC® Plus revision to ISO 6474-2; Sleeves are made of Ti-6Al-4V according to ISO 5832-3)

Inlays made of Biolox® delta ceramics may only be combined with Biolox® Delta femoral heads.

Inlays made of ELEC®Plus ceramics may only be combined with ELEC®Plus femoral heads.

Bipolar heads should be used only following a strict indication by the physician.

Combinable Ball Heads

	Article Description	Ø-Ballhead	Neck Length
	BIOLOX® Delta ISO 6474-2	28	S / M / L
		32	S / M / L / XL
		36	S / M / L / XL
		40	S / M / L / XL
	BIOLOX® Option ISO 6474-2 / ISO 5832-3	28	S / M / L / XL
		32	S / M / L / XL
		36	S / M / L / XL
		40	S / M / L / XL
	ELEC® ISO 6474-1	28	S / M / L
		32	S / M / L
		36	S / M / L
	ELEC® plus ISO 6474-2	28	S / M / L
		32	S / M / L / XL
		36	S / M / L / XL
	ELEC® plus Revision ISO 6474-2 / ISO 5832-3	28	S / M / L / XL
		32	S / M / L / XL
		36	S / M / L / XL
	Metal Ballhead CoCrMo ISO 3832-4	28	S / M / L / XL
		32	S / M / L / XL
	Metal Ballhead Steel ISO 3832-1	28	S / M / L / XL
		32	S / M / L / XL

The instrument set supplied by ARTIQO is designed for use with the implants; other instruments must not be used. (The general surgical instruments are an exception to this rule).

4.1 Application Instructions for the implantation of the hip stems

In due consideration the biomechanical principles and preoperative planning, for cementless stems the largest possible stem is implanted using the press-fit technique. The final rasp size corresponds to the stem size of the implant.

Instructions for cemented implantation:

Careful cleaning / degreasing of the medullary cavity is important for cemented application. Before cement is injected, the medullary cavity should be cleaned with a Jet-Lavage or similar and closed with an autogenous bone block or an medullary block (approx. 10-15 mm below the tip of the stem). It is recommended to use a highly viscous bone cement (PMMA), which under vacuum is mixed and inserted (3 generation cementing technique). For the centering of the stem it is advantageous if the insertion into an already viscously cement mass (at high viscosity quality approx. 3 minutes after mixing).

In the undersize-technique, to achieve the 2 mm cement mantle thickness, the stem size of the implant is selected one size smaller than the final rasp size.

The surgical technique, instructions for use and generally applicable guidelines for the implantation of cemented primary hip stems must be observed.

4.2 Application instructions for the assembly of the femoral head / Inlays of the acetabulum

The femoral stem cone and the inner cone of the femoral head have to be clean, dry and intact during the fitting. Prior to attaching the final femoral head, the cone must be thoroughly cleaned.
The appropriate femoral head is placed manually (with light axial pressure and simultaneous rotation) on the stem cone and fixed to the pole of the head with the head insertion instrument and an adequate hammer blow on the cone.

Attention! Never blow with a metal hammer on a ceramic femoral head!

When using ELEC®Plus revision ball heads, the following application notes must also be observed:

ELEC®plus revision femoral heads with titanium sleeve are supplied in separate peel pouches in one package and must be assembled before the operation. The information on the label must be observed with regard to sizes and dimensions.

ELEC®plus revision femoral heads with Ti-sleeve are intended for cup revisions, revisions after breakage of a ceramic component or when changing a ball head if a prosthesis stem left in situ, if used again a ceramic ball head. ELEC®plus revision femoral heads with Ti-sleeve can also be used for neck length corrections in primary operations.

ELEC®plus revision femoral heads with Ti-sleeve may only be used for fixation on stems that have been approved by the endoprosthesis manufacturer. Before the ELEC®plus revision femoral head can be placed on the prosthesis stem, the ceramic head must be assembled with the Ti-sleeve.

For this purpose, the Ti-sleeve is removed from the peel bag and held at its largest diameter with two fingers.

The titanium sleeve is inserted into the cone of the femoral head with a slight rotational movement until it is firmly seated. When assembling the components, make sure that they are absolutely clean and dry!

The following must be observed when placing the ELEC®plus revision femoral head on a prosthesis stem remaining in situ:

- The removal of the femoral head to be replaced must be carried out with a suitable extraction instrument to avoid damage to the stem cone remaining in situ.
- Thoroughly rinse the shaft cone remaining in situ with water to remove tissue, bone fragments or cement residues.
- Inspection of the shaft cone remaining in situ for damage and, if necessary, assessment of the damage according to the following criteria:
 - Permissible damages for the use of an ELEC®plus revision femoral head: Fine traces of damage on the used stem cone due to the removal of the previous femoral head, minor scratches up to 0.2 mm width and depth.
 - Impermissible damages for the use of an ELEC®plus revision femoral head: Crushed cone areas, flattening on the cone surface or deformations, visible notches, bevelled cone.
- Fit the ball head according to the general instructions for use of ceramic ball heads.

When using Biolox® Option-revision ball heads, the application instructions in the specific instructions for use must also be observed:

The inner cone of the acetabulum and of the outer cone of the inlay have to be clean and intact during the fitting. Prior to inserting the inlay, the inner cone must be thoroughly cleaned.

Indications

- Patients with advanced degeneration of the hip joint due to degenerative and post-traumatic osteoarthritis or rheumatoid arthritis
- Avascular necrosis of the femoral head
- Acute traumatic fractures of the femoral head or femoral neck as long as the stable anchorage of the implant is guaranteed

Contraindications

- Acute or chronic infections, whether local or systemic
- Severe muscle, nerve or vascular diseases endangering the extremity concerned
- Fractures that affect the base of the femoral neck or radiate beyond it into the trochanteric region or calcar
- Missing or poor bone substance that threatens the stable fit of the prosthesis
- Previous surgical procedures that no longer guarantee the intended support
- Marked coxa valga with a femoral neck angle >145°
- Marked coxa vara with a femoral neck angle <120°
- Any underlying condition that might compromise the function of the implant
- Revision with extensive bone defects

Precautions / limitations

- Dysplasia coxarthrosis
- Severe antetorsion of the femoral neck
- Wide femoral neck variants
- BMI>30
- Components made of ceramic are generally more resistant to wear and tear, but are more fragile

Possible side effects / complications

The adverse effects listed below are the most typical and most commonly observed consequences of hip replacement surgery.

- Loosening, degeneration and breakage of the implant components, especially at high loads
- Early or late stage infections that require removal of the implant.
- Pain, dislocation, partial dislocation, reduced range of motion, and / or instability, in particular in the event of incorrect placement of the implant or loosening or wear of the components.
- Shortening or lengthening of the leg
- Excessive wear of the polyethylene components due to intraoperative damage of the femoral spherical head, loose cement and/or bone fragments and/or an excessively high activity level or body weight of the patient.
- Femoral fractures
- Reduced joint stress
- Heterotopic ossification
- Bursitis
- Implant corrosion
- Metallosis
- Implant migration, with or without loosening of the implant
- Development of osteolysis
- Impingement
- Allergies to material components
- Venous thrombosis and pulmonary embolism
- Cardiovascular and pulmonary disorders (e. g. fat embolism)
- Haematomas
- Nerve damage

All of these complications should be diagnosed and treated as early as possible and may require revision surgery.

Risk factors for the side effects and complications listed above arise due to patient individual factors such as activity level, weight, neuromuscular status, comorbidities, traumatic events and surgical factors such as size selection, component positioning and/or axis alignment, cementing technique.

The following conditions can impair the surgical result:

- Severe osteoporosis or osteomalacia
- Severe malformations, congenial hip dislocations
- Local bone tumours
- Systemic diseases and metabolic disorders
- Falls
- History of infection
- Drug dependence or drug abuse
- Obesity
- Strenuous physical activities and those associated with strong vibrations in which the prosthesis is subjected to blows and/or excessive strain (e. g. heavy physical work, long-distance run etc.)

Prevention/minimisation of side effects

Important:

When bone cement is used, the instructions for use provided by the cement manufacturer must be followed. The implant bed must be rinsed prior to inserting the cement (for cemented fixation). Care must be taken to remove all loose particles (e. g. bone splinters) from the prepared implant bed. Following the cementing step, all excess or loose cement particles must be removed from the wound area.

If the implantation of this hip endoprosthesis is considered to be the best solution for the patient and if some of the circumstances mentioned above apply to the patient it is particularly important to brief the patient about the expected implications of these circumstances on the success of the surgical procedure.

It is also recommended to inform the patient of which activities he can undertake in order to reduce the effects of these aggravating circumstances. All information provided to the patient should be documented in writing by the surgeon performing the operation.

The surgeon must be well versed in the fixation method of the implants and the surgical technique.

Selecting the correct implant size and the correct position of the implant lowers the probability of early postoperative migration.

The use of antibiotic-substituted-bone-cement may facilitate a reduction of periprothetic infections following primary joint-replacement^[1].

Postoperative treatment

Recognised procedures should be followed for postoperative care and treatment. Documentation of the postoperative treatment should be done in accordance with internal hospital instructions and regulations, for example:

- Surgical protocol
- Postoperative X-ray images
- Patient ID card
- Regular check-ups/follow-ups

Interactions

The A2® femoral stem shaft, the ceramics and the metal heads are paramagnetic and thus react only weakly to strong magnetic fields^[2]. The higher the field strengths, the stronger the interactions. Artifacts are still to be expected in MRI examinations. However, the products were not explicitly tested according to ASTM F2503.

Patients must be instructed to avoid exceedingly strong magnetic fields and told that artificial hip joints may trigger metal detectors (e. g. at airport checkpoints).

5. Handling, maintenance and reprocessing of instruments

Section 5 applies only for instruments - Implants must never be reprocessed.

The recommendations below are for information purposes only. The manufacturer cannot be held liable for any cleaning and sterilisation processes of instruments that are performed within the purchaser's facilities. The instrument set is supplied "non-sterile" and must be cleaned, disinfected and sterilised with validated procedures prior to use.

Brief instructions for reprocessing instruments:

Our instructions are based on validation processes that were performed within the scope of reprocessing the medical device for reuse. The person concerned with reprocessing is responsible for ensuring that the actual reprocessing performed with the equipment and personnel available at the reprocessing facility achieves the desired results.

Pre-cleaning:

Prior to every automated cleaning, pre-cleaning should be carried out as soon as possible after use, preferably in an ultrasonic bath. Some washer-disinfectors are unable to reprocess instruments residue-free, especially if surgical residues have been left to dry for prolonged periods.

Automated cleaning and disinfection

The subsequent automated cleaning and disinfection should be performed with a procedure validated in accordance with the applicable rules in a cleaning and disinfection device according to EN ISO 15883 -1 and -2, preferably using (mild) alkaline cleaning agents.

Sterilisation

Steam sterilisation with a procedure validated according to EN ISO 17665 in a steam steriliser with a fractionated vacuum process according to the requirements of EN 285. The pre-validation was performed with the following program specifications:

Cycle type	Temperature	Pressure	Sterilisation time	Drying time
Fractionated vacuum	134 °C *	3050 mbar abs.	5 minutes	25 minutes

* Sterilisation temperature of 134°C (273°F) plus tolerance according to ISO 17665.

Please note the nationally applicable instructions.

Other sterilisation methods and cycles may also be used; however, these must be validated prior to implementation in routine operation according to applicable rules.

Inspection, maintenance and testing:

Dry wet or moist devices.

Repeat the cleaning and disinfection of products with residual contamination.

After each cleaning, disinfection and drying cycle, check the instruments for dryness, cleanliness and damage.
Some instruments are equipped with colour marks made from ceramic which gradually degrade with each reprocessing cycle. Check the devices for visible or palpable changes of the colour marks (cracks, detachment, loosening, protrusions or other).

 If the device no longer functions as it should, or if wear and tear of an instrument or a change of the colour mark (cracked, loose, slightly protruding colour marks) is visible, the affected parts must be either serviced or replaced.

Instruments that can be taken apart should be assembled. Only instruments that have been inspected, deemed functional and are ready for use should be placed in storage.
Ensure compatibility with the respective devices (couplings).

Packaging and transport:

Please use our transport packaging for any returns and make sure to enclose our cleaning and sterilisation form.

For more detailed instructions concerning the preparation, pre-cleaning, disinfection and sterilisation / re-sterilisation, inspection & maintenance, packaging and storage of our instruments, please refer to our reprocessing instructions (NA30-04-05).

6. Further information

For more information on the use of this product, please contact your local sales office.

7. Explanation of the Symbols Used

	Do not reuse
	Do not reprocess / do not re-sterilise
	Non-sterile
	Gamma-sterilised
	Radiation sterilised indicator (changes from yellow to red/violet for gamma treatment)
	Batch number
	Serial number
	Manufacturer
	Do not use if primary packaging is damaged
	Protect from sunlight
	Protect from moisture
	CE Mark (if applicable with 4-digit identification number of the Notified Body)
	Temperature limitation Storage temperature
	Observe warnings and instructions for use
	Observe instructions for use
	Expiry date
	Order number
	Quantity

[1] Sebastian et al. 2020: Antibiotic containing bone cement in prevention of hip and knee prosthetic joint infections: A systematic review and meta-analysis. J Orthop Translat. 2020 May 8;23:53-60. doi: 10.1016/j.jot.2020.04.005.

[2] Frank G Shellok: reference manual for magnetic resonance safety, implants and devices. (Yearly updated). Detailinformationen: <http://www.mrisafety.com/>