1. Allgemeine Informationen

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^[1] Frank G Shellock: reference manual for magnetic resonance safety, implants and devices. (Yearly updated). Detailinformationen: <http://www.mrisafety.com/>

1. General Information

ELEC[®]plus ceramic head and ELEC[®]plus revision ceramic head with titanium sleeve

Before using these products, the user is obliged to carefully study the following recommendations and notes as well as the product-specific instructions. The person placing these products on the market does not accept any liability for direct or consequential damage resulting from improper use or handling, in particular from noncompliance to the following instructions for use or from improper care or maintenance. These implants may only be used by surgeons with appropriate experience and practice in hip endoprosthetics. This product serves as a component of an artificial hip joint and cannot be used alone, but only in combination with other suitable system components.

Product description and package contents

- The product, the package contents and the materials used are specified on the product label. This product must be used in combination with the respective product-specific instructions for use and the surgical technique of the approved system components.
- The batch and SN number(s) of the implants used must be documented in the patient's records. The corresponding labels for this purpose are enclosed in the packaging of the sterile implants.

Instructions on reuse

- It is not permitted to reuse an implant which has previously been implanted in the body.
- Damage such as scratches, nicks or impact marks can lead to breakage.
- The implants are supplied sterile (sterilised by gamma radiation). They must not be re-sterilised.

Warnings / precautions

- Warnings on the packaging must be strictly observed.
- This product may only be used in its original condition. Modifications of any nature as well as mechanical manipulations are not permissible.

2. Storage and Handling

Storage of the sterile implants

Implants must always be stored in their unopened original packaging. The implants should be stored in a dry and clean place, protected from direct sunlight at a storage temperature between 10° C and 30° C. The materials were selected such that short-term temperature fluctuations, which can occur during transportation, do not have a negative impact on the performance of the implants and their packaging.

Handling of the sterile implants

⚠ Prior to insertion of the sterile implant, the packaging should be inspected for damage. If the primary packaging is damaged, the sterility is no longer guaranteed and the product must not be used.

The implants are vacuum-packed in a triple peel-bag packaging (sterilization barrier system). If it is recognizable that air has penetrated into the sterilization barrier system (air-drawn air), the sterility is also no longer guaranteed – the product is not allowed to use.

When the implant is unwrapped, a check is to be carried out to ensure that it corresponds to the designation on the packaging (article no. and size). The appropriate aseptic regulations for surgical personnel must be observed during removal of the implant from the packaging: The outermost and, if applicable, the middle PE bag are opened by the non-aseptic surgical assistant; the inner bag is opened by the aseptic surgical assistant. The implant must not come into contact with any objects which could damage its surface. Prior to insertion, every implant must be visually inspected for damage.

The ELEC[®]plus revision femoral head is supplied with a separately packed titanium sleeve (separate peel pouches in one package). The revision femoral head and titanium sleeve must be combined before fitting them onto the stem cone.

3. Instructions for Use – Special Application

Instructions Information of the patient

Prior to the surgical procedure, each patient must be informed in detail about the procedure and about possible risks.

A patient receiving a hip replacement should be informed that the life of the implant will depend on the patient's weight and how active the patient is in everyday life.

Preoperative planning

The preoperative planning provides important information about the suitable type of component, placement and possible component combinations. Additional implants should be kept handy in case other sizes are required or the intended implant cannot be used.

Restriction regarding combinations

⚠ Only implants that are supplied by ARTIQO GmbH and which have been approved for combined use can be combined with one another for hip replacement procedures. Unsuitable material or product pairings can lead to breakage, premature wear, loosening of the prosthesis, contact corrosion and other consequences. The surgeon must always ensure that the individual implant components are compatible with one another and observe the general restrictions concerning the combination.

Only the instruments of the approved system components are to be used for implantation. (The general surgical instruments are an exception to this rule). Before using the instruments, the respective enclosed instructions for use must be acknowledged. The cone dimensions of the ELEC[®]plus implants must match those of the combination implant.

ELEC[®]plus and ELEC[®]plus revision ceramic head (ISO 6474) S/M/L/XL from ARTIQO GmbH are fitted with a 12/14 cone. The cone size is marked on the label and, if possible, also on the implant. Only brand-new implants may be used.

Used materials:

The ceramic heads consist of mixed ceramic Al₂O₃/ZrO₃ (ELEC[®]plus and ELEC[®]plus revision to ISO 6474-2); Sleeves are made of Ti-6Al-4V (according to ISO 5832-3).

	ELEC [®] plus ISO 6474-2	28	S/M/L
		32	S/M/L/XL
		36	S/M/L/XL
	ELEC [®] plus Revision ISO 6474-2 / ISO 5832-3	28	S/M/L/XL
		32	S/M/L/XL
		36	S/M/L/XL

Surgical technique for the assembly of femoral heads:

In order to ensure a safe function of the implant, the utmost care taken during implantation is an indispensable requirement. For primary treatments, remove the protective cap, which protects the stem cone from damage, just before adjusting the femoral head. The following must be observed before placing the femoral head onto the stem cone:

Thoroughly rinse the cone with water to remove any remaining tissue, bone fragments or cement residues.

The femoral stem cone and the inner cone of the femoral head have to be clean, dry and intact during the fitting. Prior to attaching the final femoral head, the cone must be thoroughly cleaned. The appropriate femoral head is placed manually (with light axial pressure and simultaneous rotation) on the stem cone and fixed to the pole of the head with the head insertion instrument and an adequate hammer blow on the cone.

- Never blow with a metal hammer on a ceramic femoral head!
- Once an ELEC[®]plus femoral head has been fitted on a stem cone and removed again, it must not be reapplied again. Instead, a metal head or an ELEC[®]plus revision femoral head must be used. The same applies if the prosthesis stem is left in situ during revision but the femoral head is replaced.

When using ELEC[®]plus revision ball heads, the following application notes must also be observed:

ELEC[®]plus revision femoral heads with titanium sleeve are supplied in separate peel pouches in one package and must be assembled before the operation. The information on the label must be observed with regard to sizes and dimensions

ELEC[®]plus revision femoral heads with Ti-sleeve are intended for cup revisions, revisions after breakage of a ceramic component or when changing a ball head if a prosthesis stem left in situ, if used again a ceramic ball head.

ELEC[®]plus revision femoral heads with Ti-sleeve can also be used for neck length corrections in primary operations.

ELEC[®]plus revision femoral heads with Ti-sleeve may only be used for fixation on stems that have been approved by the endoprosthesis manufacturer

Before the ELEC[®]plus revision femoral head can be placed on the prosthesis stem, the ceramic head must be assembled with the Ti-sleeve. For this purpose, the Ti-sleeve is removed from the peel bag and held at its largest diameter with two fingers. The titanium sleeve is inserted into the cone of the femoral head with a slight rotational movement until it is firmly seated. When assembling the components, make sure that they are absolutely clean and dry!

The following must be observed when placing the ELEC[®]plus revision femoral head on a prosthesis stem remaining in situ:

- The removal of the femoral head to be replaced must be carried out with a suitable extraction instrument to avoid damage to the stem cone remaining in situ.
- Thoroughly rinse the shaft cone remaining in situ with water to remove tissue, bone fragments or cement residues.
- Inspection of the shaft cone remaining in situ for damage and, if necessary, assessment of the damage according to the following criteria:
 - Permissible** damages for the use of an ELEC[®]plus revision femoral head: Fine traces of damage on the used stem cone due to the removal of the previous femoral head, minor scratches up to 0.2 mm width and depth.
 - Impermissible** damages for the use of an ELEC[®]plus revision femoral head: Crushed cone areas, flattening on the cone surface or deformations, visible notches, bevelled cone.
- Fit the ball head according to the general instructions for use of ceramic ball heads

Indications

- Total endoprosthesis, primary surgery in combination with prosthetic stems and acetabular cups
- Hemiarthroplasty
- Revision surgery (= endoprosthesis removal) using new prosthesis stems and cups
- Partial revision if a prosthesis stem left in situ without damage or with permissible damage, only permissible with ELEC[®]plus revision ceramic head

Contraindications

- Revision leaving shafts with impermissibly damaged shaft cones in place (danger of component fracture)
- Use with non-compatible system components

Possible side effects / complications

The adverse effects listed below are the most typical and most commonly observed consequences of hip replacement surgery.

- Dislocation, subluxation, reduced range of motion, unwanted change in leg length .
- Infections, haematomas, delayed wound healing
- Venous thrombosis and pulmonary embolism and other cardiovascular disorders.
- Displacement or loosening of the implant
- Fracture of implant, bone or cement
- Side effects due to sliding bearing wear (e.g. osteolysis, inflammatory tissue reactions)

The side effects and undesirable effects listed below can particularly occur in connection with ceramic components such as the ELEC[®]plus and ELEC[®]plus revision femoral head:

- Material failure, ceramic fracture
- Noise generation (e.g. squeaking, clicking)
- Corrosion at the cone connection

The femoral head serves as a component of an artificial hip joint. Therefore, the possible side effects / complications in the instructions for use of the other system components must also be taken into consideration.

Side effects should generally be diagnosed and treated as early as possible.

Prevention / minimisation of side effects

Important:

If the implantation of this hip endoprosthesis is considered to be the best solution for the patient and if some of the circumstances mentioned above apply to the patient it is particularly important to brief the patient about the expected implications of these circumstances on the success of the surgical intervention.

Patients receiving a hip replacement must be advised that the life of the implant depends on their weight and activity level.

It is also recommended to inform the patient of which activities he can undertake in order to reduce the effects of these aggravating circumstances. All information provided to the patient should be documented in writing by the surgeon performing the operation.

These implants may only be used by surgeons with appropriate experience and practice in hip endoprosthetics. The surgical technique must be known to the surgeon.

Postoperative treatment

Recognised procedures should be followed for postoperative care and treatment. Documentation of the postoperative treatment should be done in accordance with internal hospital instructions and regulations, for example:

- Surgical protocol
- Postoperative X-ray images
- Patient ID card
- Regular check-ups/follow-ups

Interactions








Ceramic ball heads and the titanium sleeve are paramagnetic and thus react only weakly to strong magnetic fields^[1]. The higher the field strengths, the stronger the interactions. Artifacts are still to be expected in MRI examinations. However, the products were not explicitly tested according to ASTM F2503.

Patients must be instructed to avoid exceedingly strong magnetic fields and told that artificial hip joints may trigger metal detectors (e. g. at airport checkpoints).

4. Further Information

For more information on the use of this product, please contact your local sales office.

5. Explanation of the Symbols Used

	do not reuse
	do not reprocess / do not re-sterilise
	Gamma-sterilised
	Batch number
	Serial number
	Manufacturer
	do not use if primary packaging is damaged
	protect from sunlight
	protect from moisture
	CE Mark according to Directive 93/42/EEA
	Temperature limitation / Storage temperature
	Observe warnings and instructions for use
	Observe instructions for use
	Expiry date
	Order number
	Quantity
	Radiation sterilised indicator (changes from yellow to red/violet for gamma treatment)

^[1] Frank G Shellock: reference manual for magnetic resonance safety, implants and devices. (Yearly updated). Detailinformationen: <http://www.mrisafety.com/>