OSARTIS®

POLY PLUG 18 / POLY PLUG 24

EN

POLY PLUG 18 / POLY PLUG 24 INSTRUCTIONS FOR USE

Screw on the intramedullary restrictor on the applicator with a clockwise rotary motion. Warning: be sure not to cross-thread or over-tighten the end cap of the applicator. Insert the restrictor into the medullary canal and position it approximately 2 centimetres lower than the tip of the stem. Once inserted, disengage the inserter by twisting counter-clockwise. If a secondary brush and lavage is desired, take care not to dislodge the cement restrictor.

WARNINGS

Strictly aseptic surgical techniques must be observed. Sterility is assured only if the unit container is not damaged or opened. The device is contraindicated for non-cemented total hip procedures and for applications other than those which are indicated. Moreover it is contraindicated in infectious arthritis and in active infection of the joint(s) to be replaced. Use of the cement restrictors is also contraindicated where the loss of musculature or neuromuscular compromise in the affected limb would render the surgical procedure unjustifiable. The insertion of the restrictor should be done with caution to avoid the risk of pulmonary embolism (BCIS).

CAUTION

The device is single-use and intended for single patient. Do not re-sterilize and/or re-use any of the components. The device may not be reused as there is a risk of infection for the patient and the same performance cannot be guaranteed. In addition, resterilization may alter the morphology and/or the chemical/physical/mechanical properties of the device eventually causing its malfunctioning with serious risks for the patient's health. All residues must be considered waste surgical material and therefore should be eliminated at the end of the operation.

STORAGE AND TRANSPORT CONDITIONS

- Store and transport at a temperature below 25 °C and in a cool and dry environment. However, if the temperature during the transport cannot be maintained under 25 °C with a relative humidity < 90%, the medical device can tolerate a temperature between 0 °C and 40 °C with a relative humidity less than 90% for maximum 5 days.
- The device must be shipped and transported in a double wall corrugated box with ECT $\geq$ 44 lbf/in (or ECT $\geq$ 7.71 kN/m).

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Osartis GmbH
TECRES S.p.A.