

Surgical technique

4-MOTION®

Total knee replacement system



4-motion®

4-motion® TOTAL KNEE SYSTEM


Table of contents

4-motion® knee endoprosthesis	4-5
Design features	6-7
Femoral component	6
Tibial component PSA	7
Indications and contraindications	8
Customised instrument concept	9
Workflow / Process	10
Preoperative planning	11-12
4-motion® online planning tool	11
Recommendations for the positioning of the implants and determination of the size	12
Implantation	13-20
Introduction and preparation	13
Distal femoral resection	13
Implantation	14
Femoral AP and chamfer resection	15
Final femoral preparation	16
Tibial resection	17
Tibial preparation	18-20
Trial repositioning	20
Ligament-balancing technique	21-24
Description of the function of the 4-motion® ligament balancer	21
Ligament-guided surgical technique	22-24
Implantation of the final implants	25-27
Cemented / non-cemented anchoring of the tibial and femoral components	25
Implantation of the patellar component	26
Compatibility table	26
Insert versions	27
Postoperative treatment / Sterilisation	28
Implants / Range	29-35
Instruments / Range	36-37



4-motion®

K N E E E N D O P R O S T H E S I S



With the development of the 4-motion® knee endoprosthesis, we are consistently pursuing the reconstruction of the anatomical planes of movement in the knee joint. This reconstruction offers numerous advantages in the context of prosthesis implantation. The femoral offset can be reconstructed since the operation-related external rotation of the femoral component is no longer necessary. In addition, the anatomical position significantly reduces the need for ligament releases.

Design-Merkmale



Femoral component

- Anatomical joint line
- Reconstruction of the dorsal offset
- Q angle of 9°
- Widened patella inlet
- Different dorsal condyle length
- 10 sizes in 3 mm a-p steps
- 2 different m-l dimensions for the relevant sizes





Tibial component PSA

- Asymmetric base plate
- Central keel with semicircular shield
- Reconstruction of the patient-specific anatomical alignment in the range of 87° to 90° (MPTA) possible
- Suitable for allergy-sensitive patients (titanium alloy)



Indications and contraindications

Indications

- Patients with advanced degeneration of the knee joint resulting from degenerative and post-traumatic arthritis (PTA) or rheumatoid arthritis (RA)
- Avascular necrosis of the femoral condyle
- Moderate varus, valgus, or flexion deformity

Contraindications

- Acute or chronic infections, local or systemic
- Severe muscular, neural, or vascular diseases putting the affected extremity at risk
- Absent bone substance or inadequate bone quality posing a risk to the stable fit of the prosthesis
- Severe osteoporosis or osteomalacia
- Hypersensitivity to the materials employed
- Loss of the ligaments
- Any concomitant disease capable of putting the function of the implant at risk
- Revision with pronounced bony defects

Precautions / restrictions

- Major malpositions in combination with instabilities
- BMI >30
- Loss of posterior cruciate ligament
- High level of physical activity (e.g., competitive sports, strenuous physical work)
- Diabetes

Customised instrument concept

“Image to instrument” procedure


This innovative procedure allows for customized endoprosthetic knee replacement and begins with a CT scan of the patient, which is uploaded to a secure, web-based planning platform operated by (<https://artiqo.mediation-planning.com>).

On the basis of these DICOM data, an interactive, three-dimensional CAD model of the patient's pathologically altered knee is generated, and a detailed, preoperative proposal sent to the surgeon for review and approval. The important aspect here is that the surgeon always has the final decision regarding the procedure. Depending on the anatomical requirements, the user can adjust the axial and rotational alignment parameters as required. Once the design has been finalised, patient-specific cutting blocks are generated. These enable the surgeon to implement the highly precise, three-dimensional planning during the operation.


This “image to instrument” concept combines various functions and offers many advantages for both surgeons and patients alike:

- 3D reconstruction of the knee to be operated on
- Interactive 3D web planning
- Precise positioning of the implant
- No need to open the medullary canal
- Fewer surgical steps and shorter operation time
- Reduction of time required and costs for cleaning, assembly, and sterilisation

The instruments manufactured on this basis are precisely adapted to the individual surface contour and anatomy of the respective knee joint.




Workflow and process



MEDivation is pleased to be the contact partner for the commissioned radiology centres and provides them with the technical information for the generation of the patient's image data, so that the scans can always be performed in accordance with the required protocol and that at the first attempt.

In principle, CT scans are planned in three subareas (hip, knee, and ankle joint). For better orientation and increased precision, a whole leg radiograph can also be included in the planning. Once the data are available, MEDivation prepares the segmentation and planning proposal. On the basis of the images, the resection planes and the positioning of the components are precisely aligned with the mechanical leg axes and the clearly identified landmarks. Osteophytes are taken into consideration during the design of the cutting blocks, as they are an important part of the current patient anatomy and essential for the precise fit of the instrument to the bone. A lead time of approx. seven weeks should be planned for the entire cycle.



Preoperative planning



4-motion® online planning tool


The preoperative planning begins with the input of the patient data, specification of the desired time and date of the surgery and sending of the CT data to MEDivation's web-based planning platform (<https://artiqo.medivation-planning.com>).

NOTE:

A surgeon-specific procedure such as the femur-first or ligament-supported tibia-first technique can be predefined directly in the customer profile.

The automated planning utilises a planning algorithm which takes the data regarding the leg axis (varus, neutral or valgus), pathology (tibia and / or femur) and degree of instability into consideration during positioning of the implants.

- Reconstruction of the joint line between 87° and 90° depending on the anatomy
- For varus deformities and neutral leg axes, the preoperative alignment (LDFA) of the femur is reconstructed
- For valgus deformities, the preoperative alignment (MPTA) of the tibia is reconstructed
- 4° slope
- Neutral flexion position of the femur in the sagittal plane
- Rotational orientation of the component on the posterior femoral condyle
- Determination of the size and a-p position through primary orientation to the ventral cortex of the femur



NOTE:

Depending on the anatomical requirements, the surgeon can check and individually adjust the axial and rotational orientation parameters in the planning proposal. Further details on the workflow can be found in MEDivation's "Webplanning" user manual.

Finally, a summary of the implantation parameters and geometries of the patient-specific cutting blocks is prepared for approval.

Preoperative planning



Frontal plane

Recommendations for the positioning of the implants and determination of the

Femoral component

The femoral component must sit flush against the cortex without protruding over the bone beneath. In case of doubt, the mediolaterally smaller implant is selected.

In the frontal plane, the components are aligned in the anatomical joint line. Alignment in the range LDFA $87^\circ \pm 3^\circ$ (depending on the anatomy and defect situation) is recommended.

Protrusion of the patellar sliding surface in the ventral direction must always be avoided in order to minimise retropatellar contact forces and achieve good flexion.

The rotation is aligned parallel to the posterior femoral condyles as standard.


Tibial component

The tibial component must sit flush against the cortex without protruding over the bone beneath.

In the frontal plane, the components are aligned in the anatomical joint line. Alignment in the following range (depending on the anatomy, tibial implant type and defect situation) is recommended:

■ Tibia PSA: MPTA 87° to 90°

Special attention should be paid to the correct rotational position of the tibial component. Alignment over the attachment point of the posterior cruciate ligament to the medial third of the tibial tuberosity is recommended.



Horizontal plane

Implantation


Introduction and preparation

The surgical steps described below are based on the femur-first technique, which is the standard option on the planning platform for the **implantation of the 4-motion® knee endoprosthesis**.

The 4-motion® knee endoprosthesis can optionally also be implanted using a tibia-first or ligament-supported tibia-first technique depending on the preferred surgical procedure. **Detailed information can be found in the “Ligament-balancing technique” section.** Instruments for the ligament-guided surgical technique are not included in the instrument tray as standard. If required, the customer profile on the planning platform also offers the option of individually predefining the operator-specific procedure or specifying the surgical technique or requirement when placing the order.

NOTE:

The patient-specific cutting blocks are custom-made devices produced to fit the patient's anatomy based on the operating surgeon's preoperative planning. The user must ensure that the correct cutting blocks for the patient undergoing surgery are used. The patient code (initials, date of birth, side) is printed on the packaging and on all instruments. The cutting blocks are designed for a 1.27 mm saw blade and pins with a diameter of 3.2 mm as standard and may only be used with the specified implant system. Use of other instruments can result in damage to the cutting blocks and undesirable outcomes.





Distal femoral resection

The patient-specific femoral cutting block is positioned stably on the bone at the correct contact sites as per the planning diagram. The edge of the articular surface with the osteophytes must be clear. Soft tissue and cartilage at the contact points should be removed in order to be able to position the cutting block precisely.

NOTE:

Special attention must be given to ensuring that all cartilage is removed from the contact area for the cutting block. It is recommended to position the cutting block temporarily, mark the contact area and then remove the cartilage with a curette. Osteophytes should not be removed unless their presence would prevent the correct positioning of the block. Sterilizable plastic phantoms of the femur can be optionally ordered for trial positioning, enabling identification of the correct position.




Implantation

The femoral cutting block is fixed in position from the ventral direction with drill pins. The subsequent fixation holes for the 4in1 cutting block are then predrilled at least 30 mm deep via the distal guide holes using a 3.2 mm drill bit.

NOTE:


Make sure the cutting block does not slip and is not displaced during the fixation. Verify the correct positioning again once the pins are in place.



The distal part of the cutting block must be removed before the bone is cut.

NOTE:


Ensure that the distal drilling holes for the 4in1 cutting block have been drilled before removal.



The distal cut can then be performed directly through the cutting block.

NOTE:

If necessary, a correction cutting block can be used for parallel recutting or a varus / valgus correction of 2° performed.



Femoral AP and chamfer resection

Position the 4in1 cutting block on the inserted 3.2 mm pins. When doing so, ensure that the cutting block sits level and flush against the already resected distal surface. The final fixation is performed with 3.2 mm oblique pins.

NOTE:

To avoid errors when positioning the pins, it is recommended to insert the pins into the predrilled holes with the cutting block level in place.

Ensure that the planned size of cutting block is used. (Cutting blocks in the next sizes up and down can also be optionally ordered).

Check the cut planes with the cortical palp and make the anterior and posterior bone cuts. The distal pins must be removed for the chamfer cuts.

NOTE:

If required, the cutting block can be moved ± 2 mm in the anterior/posterior direction. For precise bone cuts, ensure that the saw blade is inserted smoothly through the saw slot with no lever action.



Final femoral preparation


The correct femoral trial component is positioned on the femur, using the femoral impactor if necessary.

NOTE:

The femoral component must sit flush against the cortex without protruding over the bone beneath. The general recommendation is to align the femoral component in a lateral direction.

The final positioning in the m/l direction is determined based on the patellar tracking.




The definitive m/l position of the femoral component is fixed via the drilling of the plug holes. The central axis for the trochleaplasty is then marked with a pen or electrocautery.



Finally, the trochleaplasty is performed in a V shape using the rasp. First, the side of the rasp with the narrow ridge is placed on the marked line and the ridge rasped as far as it will go. The rasp is then turned over and the full toothed surface of the broad side applied along the pre-rasped surface.

NOTE:

Ensure that the trochea rasp is guided in the plane of the ventral chamfer.



Tibial resection

The patient-specific tibial cutting block is positioned stably on the bone at the correct contact sites as per the planning diagram.

NOTE:


Special attention must be given to ensuring that all cartilage is removed from the contact area for the cutting block. It is recommended to position the cutting block temporarily, mark the contact area and then remove the cartilage with a curette. Sterilizable plastic phantoms of the tibia can be optionally ordered for trial positioning.

The tibial cutting block is fixed in position with 2 parallel drill pins and an oblique pin. The cut plane is then checked with the cortical palp, the longitudinal axis verified with the alignment rod and the tibial cut performed. Ensure that the posterior cruciate ligament is preserved.

NOTE:

The tibial resection is performed with a 1.27 mm saw blade as standard. For a precise bone cut, ensure that the saw blade is inserted smoothly through the saw slot with no lever action.

If necessary, a correction cutting block can be used for parallel recutting or a varus / valgus correction of 2° performed.




Tibial preparation

The planned size is checked again in situ using the tibial base plate. Complete cortical coverage without protrusion of the components is the aim. Note that the plate can be used side-specifically from both sides for a right and a left knee.

NOTE:

The tibial size must not be more than two sizes different (larger or smaller) from that of the femur. The corresponding base plate must be used depending on the planning and selected implant type.



pecial attention should be paid to the correct rotational position. This is done using the femoral trial, tibial base plate, and trial insert. When the knee joint is actively flexed and extended, the implant finds its individually desired position. This rotational position is marked on the transition from the front edge of the tibia to the tibial base plate with an electrocautery or similar, so that it can be found again later.

Femoral trial		Marking	2s	3s	4s	5s 5w	6s 6w	7w	8w	9w
Tibial base plate	1									
	2									
	3									
	4									
	5									
	6									
	7									
	8									
	9									
	10									

NOTE:

Alternatively, the rotational alignment can also be performed anatomically over the attachment point of the posterior cruciate ligament to the medial third of the tibial tuberosity.


NOTE:

In principle, there are three trial insert sizes (A, B and C) available in four heights (6, 8, 10, 12). The compatibility of the insert groups with the femur sizes is shown in the table. In addition, the femoral and tibial base plates are marked with the colour code of the corresponding insert group.

Information on the insert range

The 4-motion® product range offers the following CR insert:

CR insert VEX for increased femoral range of motion (deep flexion and rotation) Detailed information on the insert versions can be found in the description on page 27. If required, the availability of both insert versions must be checked prior to the implantation.




Tibial preparation

Finally, the defined rotational position is checked with the alignment rod and the tibial base plate fixed in place with two pins.

NOTE:


The corresponding angle socket (0° – 3°) based on the planned MPTA angle on the tibia (90° – 87°) is inserted into the handle. When doing so, it is important to ensure that the angle socket is employed for the correct side (left and right) and that the arrowhead points towards the head of the tibia when correctly aligned.



The preparation for the tibial keel begins with the application of the drill and keel tower.

NOTE:

The corresponding keel tower (Tibia PSA) must be used depending on the planning and selected implant type.



The tibial keel is then reamed with the stepped drill.

NOTE:

There is one drill available for all sizes, which is equipped with an AO quick coupling as standard. The corresponding drill must be used depending on the planning and selected implant.

The drill is inserted as far as it will go.




Tibial preparation

The corresponding keel punch is selected based on the size of the tibial implant and connected with the prime handle. The keel punch is impacted as far as it will go.

NOTE:

The keel punch is available in three sizes: size 1-3, size 4-6 and size 7-10



Trial repositioning

Final trial repositioning with checking of the ligamentary stability and patellar tracking is recommended.


Ligament-balancing technique

Description of the function of the 4-motion® ligament balancer

1. ASSEMBLY

The ligament balancer is a four-part instrument comprising:

- A Balancer body
- B Balancer spacer
- C Extension attachment
- D Flexion attachment




The first step of assembly of the ligament balancer is insertion of the balancer spacer into the balancer body. The extension or flexion attachment can be selected and connected depending on the intended use (verification of the extension or flexion gap).

2. GENERAL FUNCTIONING

For proper functioning of knee joints, the gap should be as parallel and similarly sized in both extension and flexion. A variety of techniques can be employed to achieve this. In contrast to conventional techniques (e.g., the measured resection technique), the patient-specific soft-tissue tension and alignment are taken into consideration in ligament-guided techniques.

The 4-motion® ligament balancer can be used intraoperatively to verify the extension gap in relation to the parallelism and height of the resections and to balance the flexion gap based on the tibial resection with ligament guidance.

The aim is to find an improved alignment and soft-tissue balance of the 4-motion® anatomical knee prosthesis utilising the ligament-balancing technique and the defined workflow.




Ligament-guided surgical technique

NOTE:

The prerequisite for use of the ligament balancer in extension is the (temporary or final) resection on the distal femur and the resection on the tibia as well as the removal of all osteophytes.

The ligament balancer is assembled as shown under Point 1 “Assembly” and the extension module attached.

The ligament balancer is inserted into the joint gap with the leg axis in a neutral “not hyperextended” position. When doing so, ensure that the ligament balancer is completely closed and not spread.



Then use the hexagon socket on the side to tighten the ligament balancer in the direction of the arrow using the torque wrench provided (2 Nm).

NOTE:


The applied torque is used to limit the clamping force and for comparison of the extension and flexion gap.

The two indicators can be utilised for estimation of the parallelism and the height of the joint gap.

NOTE:

Angle differences can be resolved via a correction cut on the femur or ligament releases.

The minimum height required for the implantation of the knee prosthesis with the lowest PE height (6 mm) is 18 mm.




Ligament-guided surgical technique

Verification of the flexion gap:

NOTE:


For the use of the ligament balancer in flexion, the 4in1 cutting block is first pushed onto the pins already positioned using the PSI femur cutting block.

The ligament balancer is inserted into the joint gap at approx. 90° flexion before resection of the anterior and posterior femoral cut. When doing so, ensure that the ligament balancer is completely closed and not spread.



NOTE:


The thickness of the ligament balancer / the distance to the surface when closed is 10 mm and thus corresponds to the minimum height required for implantation of the tibial component with the lowest PE height (6 mm).



The flexion attachment is inserted into the ligament balancer and connected to the 4in1 cutting block using the two locking pins. Ensure that the pins are connected to the cutting block over at least half the pin length. If necessary, the ligament balancer must be inserted further into the flexion gap or the flexion angle of the knee joint (approx. 90°) checked and corrected.

Ligament-guided surgical technique

Then, as in extension, the ligament balancer is tightened in the direction of the arrow using the hexagon socket on the side and the torque wrench provided.



NOTE:

The locking mechanism is released (unclamping) using the button on the side of the balancer body.

Once the 4in1 cutting block has been positioned, the remaining steps are performed as described in the standard surgical technique above.

Implantation of the final implants

In principle, the 4-motion® knee system can be implanted with or without the use of cement. It is important to note that the respectively correct implants are selected for the respective technique and then the compatibility checked again.

Cemented anchoring of the tibial and femoral components:

When cementing components, the recommendation is to use two portions of cement.

The procedure usually begins with implantation of the tibial component via anterior subluxation of the tibia. The first portion of bone cement is mixed in accordance with the respective manufacturer's instructions. A layer of cement is applied to the underside of the tibial component, around the keel and to the resected tibial surface. The tibial component is inserted into the tibia by hand and impacted using the tibial impactor. Excess cement is removed thoroughly. The definitive tibial PE insert is inserted into the tibial component and fixed in place via the ventral catch using the tibial impactor. PE inserts which have already been fixed in place and removed from the tibial component again must not be reused.

For the implantation of the femoral component, the knee is flexed in an approx. 90° position and a collateral ligament retractor placed both laterally and medially. Cement is applied to the underside of the implant and the resected femoral surfaces. The femoral component is applied to the femur. When doing so, ensure that the plugs align with the holes and no soft tissues are trapped beneath the implant. The final fixation of the femoral component is performed with the femoral impactor. Scratches on the surfaces of the components must be avoided. The retractors are removed again. The medial and lateral sides are checked to ensure that the femoral component has been impacted completely. Excess cement is removed thoroughly.

Non-cemented anchoring of the tibial and femoral components:

The non-cemented procedure also begins with the tibial component. The tibia is subluxated and the tibial component inserted manually into the tibia. The final fixation is performed with the tibial impactor. The definitive tibial PE insert is inserted into the tibial component and fixed in place via the ventral catch using the tibial impactor. PE inserts which have already been fixed in place and removed from the tibial component again must not be reused.

For the implantation of the femoral component, the knee is flexed in an approx. 90° position and a collateral ligament retractor placed both laterally and medially. The femoral component is applied to the femur. When doing so, ensure that the plugs align with the holes and no soft tissues are trapped beneath the implant. The final fixation of the femoral component is performed with the femoral impactor. Scratches on the surfaces of the components must be avoided.




Implantation of the patellar component:

The underside of the patellar component is covered with cement and the three plug holes in the patella filled with cement. With the leg extended, the patellar component is inserted by hand and pressed in using the patellar forceps. Excess cement is removed. The patellar forceps are left in place until the cement has set fully.

NOTE:

Preparation of the patella (underside replacement) is performed optionally. The 4-motion® knee endoprosthesis can optionally be used with or without an artificial patella replacement. The patellar instruments are not part of this surgical technique and can be requested if required.

Compatibility overview

 Dia. 26 Dia. 28 Dia. 32 Dia. 35												
			Insert									
Patella	Femur first		1	2	3	4	5	6	7	8	9	10
	2											
	3											
	4											
	5											
	6											
	7											
	8											
	9											
			1	2	3	4	5	6	7	8	9	10
Tibia												
												

Insert versions for the 4-motion® knee endoprosthesis

The 4-motion® product range contains two CR insert versions. Both insert versions have a concave shape from anterior to posterior, which has the priority of limiting the translational movement of the femur in the anterior direction and can partially compensate for the stabilizing function of the missing anterior cruciate ligament.

The CR insert made of standard polyethylene has a raised anterior and posterior lip (compared to the PE insert VEX). This insert version can be utilised to achieve increased a-p stability for patients with lax joint conditions.

Compared to the PE insert Standard, the VEX insert has been designed for patients with good joint stability and high demands on mobility. In particular, the lower dorsal lip enables an extended range of motion for deep flexion. In addition, the smaller concavity provides increased freedom of rotation of the femur.



	CR insert Standard with raised anterior and posterior lip	CR insert VEX
Material	UHMWPE in accordance with ISO 5834-2	Highly crosslinked UHMWPE with vitamin E
Unconstrained insert type	CR insert for increased femoral a-p stability "Deep Dish"	CR insert for increased femoral range of movement (deep flexion and rotation)
Rotational freedom of the femur in extension*	6°	12°
Rotational freedom of the femur in flexion*	12°	20°
Posterior slope of the PE inserts incl. 4° slope	8,9°	8,5°
		

*Calculation for femur size 6 and insert size 4

Postoperative treatment / Sterilisation

Postoperative treatment

The postoperative treatment depends on the result of the surgery. In principle, mobilisation can be begun early. Questions such as partial or full weight-bearing, crutches, three- or four-point gait, etc., are recommended by the surgeon or hospital.

The bone quality and patient's condition must always be taken into consideration here. Physiotherapy during the stay in hospital is recommended.




Sterilisation

Implants

All implants described in the surgical technique are supplied sterile by the manufacturer. Resterilisation is not permitted.

Instruments

The system elements and instruments are not delivered sterile. They must be cleaned, disinfected and sterilized in a validated procedure before being used. Cleaning instructions contained validated instructions for the reprocessing of a medical device can be provided. The reprocessor is responsible for ensuring that the reprocessing actually performed with the equipment, materials and personnel available in the reprocessing facility achieves the desired results. Validation and routine monitoring of the process are normally required for this. Instrument manufacturers and distributors assume no liability for the sterilisation of devices by the purchaser.



Implants

4-motion® femoral components

Size left	Size right	A-P	M-L slim	M-L wide
2s-L	2s-R	53	58	–
3s-L	3s-R	56	60	–
4s-L	4s-R	59	62	–
5s-L	5s-R	62	64	–
6s-L	6s-R	65	66	–
5w-L	5w-R	62	–	68
6w-L	6w-R	65	–	71
7w-L	7w-R	68	–	74
8w-L	8w-R	71	–	77
9w-L	9w-R	74	–	80

4-motion® tibial components PSA (3° tibial keel optionally available)

Size left	Size right	M-L	A-P medial	A-P lateral
1-L	1-R	63	42	40
2-L	2-R	66	44	42
3-L	3-R	69	46	44
4-L	4-R	72	48	46
5-L	5-R	75	50	48
6-L	6-R	78	52	50
7-L	7-R	81	54	52
8-L	8-R	84	56	54
9-L	9-R	87	58	56
10-L	10-R	90	60	58

Implants

4-motion® femoral components cemented

Side	Size	Art.-no
Right	2s-R*	4801-RS02
Right	3s-R	4801-RS03
Right	4s-R	4801-RS04
Right	5s-R	4801-RS05
Right	6s-R	4801-RS06
Right	5w-R	4801-RW05
Right	6w-R	4801-RW06
Right	7w-R	4801-RW07
Right	8w-R*	4801-RW08
Right	9w-R*	4801-RW09

Left	2s-L*	4801-LS02
Left	3s-L	4801-LS03
Left	4s-L	4801-LS04
Left	5s-L	4801-LS05
Left	6s-L	4801-LS06
Left	5w-L	4801-LW05
Left	6w-L	4801-LW06
Left	7w-L	4801-LW07
Left	8w-L*	4801-LW08
Left	9w-L*	4801-LW09

* Special size only on request / optionally available

Implant material:

CoCrMo casting alloy according to ISO 5832-4

Surface: Sliding surfaces polished

NOTE:

On request, the femoral components are also available in an allergy-inhibiting “Sensitive” version with a titanium-niobium nitride (TiNbN) coating for allergy sufferers.

Implants

PSA tibial components cemented

Side	Size	Art.-no.
Right	1*	4631-R001
Right	2	4631-R002
Right	3	4631-R003
Right	4	4631-R004
Right	5	4631-R005
Right	6	4631-R006
Right	7	4631-R007
Right	8	4631-R008
Right	9	4631-R009
Right	10*	4631-R010

Left	1*	4631-L001
Left	2	4631-L002
Left	3	4631-L003
Left	4	4631-L004
Left	5	4631-L005
Left	6	4631-L006
Left	7	4631-L007
Left	8	4631-L008
Left	9	4631-L009
Left	10*	4631-L010

* Special size only on request / optionally available

Implant material of PSA tibial component:

Wrought Ti6Al4V alloy in accordance with ISO 5832-3

3° tibial components cemented (optional)

Side	Size	Art.-no.
Right	1*	4601-R001
Right	2	4601-R002
Right	3	4601-R003
Right	4	4601-R004
Right	5	4601-R005
Right	6	4601-R006
Right	7	4601-R007
Right	8	4601-R008
Right	9	4601-R009
Right	10*	4601-R010

Left	1*	4601-L001
Left	2	4601-L002
Left	3	4601-L003
Left	4	4601-L004
Left	5	4601-L005
Left	6	4601-L006
Left	7	4601-L007
Left	8	4601-L008
Left	9	4601-L009
Left	10*	4601-L010

* Special size only on request / optionally available

Implant material of 3° tibial component:

CoCrMo casting alloy in accordance with ISO 5832-4

Implants

PSA tibial components non-cemented

Side	Size	Art.-no.
Right	1*	3631-R001
Right	2*	3631-R002
Right	3*	3631-R003
Right	4*	3631-R004
Right	5*	3631-R005
Right	6*	3631-R006
Right	7*	3631-R007
Right	8*	3631-R008
Right	9*	3631-R009
Right	10*	3631-R010

Left	1*	3631-L001
Left	2*	3631-L002
Left	3*	3631-L003
Left	4*	3631-L004
Left	5*	3631-L005
Left	6*	3631-L006
Left	7*	3631-L007
Left	8*	3631-L008
Left	9*	3631-L009
Left	10*	3631-L010

* Special size only on request / optionally available

Implant material of PSA tibial component:

Wrought Ti6Al4V alloy in accordance with ISO 5832-3, TPS

Implants

4-motion® insert VEX, highly crosslinked UHMWPE with vitamin E

Art.-no.	Size	Height
7411-0106	Size 1	6 mm
7411-0107*	Size 1	7 mm
7411-0108	Size 1	8 mm
7411-0109*	Size 1	9 mm
7401-0110	Size 1	10 mm
7411-0111*	Size 1	11 mm
7411-0112*	Size 1	12 mm
7411-0115*	Size 1	15 mm
7411-0206	Size 2	6 mm
7411-0207*	Size 2	7 mm
7411-0208	Size 2	8 mm
7411-0206*	Size 2	9 mm
7411-0210	Size 2	10 mm
7411-0211*	Size 2	11 mm
7411-0212*	Size 2	12 mm
7411-0215*	Size 2	15 mm
7411-0306	Size 3	6 mm
7411-0307*	Size 3	7 mm
7411-0308	Size 3	8 mm
7411-0309*	Size 3	9 mm
7411-0310	Size 3	10 mm
7411-0311*	Size 3	11 mm
7411-0312*	Size 3	12 mm
7411-0315*	Size 3	15 mm

Art.-no.	Size	Height
7411-0406	Size 4	6 mm
7411-0407*	Size 4	7 mm
7411-0408	Size 4	8 mm
7411-0409*	Size 4	9 mm
7401-0410	Size 4	10 mm
7411-0411*	Size 4	11 mm
7411-0412*	Size 4	12 mm
7411-0415*	Size 4	15 mm
7411-0506	Size 5	6 mm
7411-0507*	Size 5	7 mm
7411-0508	Size 5	8 mm
7411-0506*	Size 5	9 mm
7411-0510	Size 5	10 mm
7411-0511*	Size 5	11 mm
7411-0512*	Size 5	12 mm
7411-0515*	Size 5	15 mm
7411-0606	Size 6	6 mm
7411-0607*	Size 6	7 mm
7411-0608	Size 6	8 mm
7411-0609*	Size 6	9 mm
7411-0610	Size 6	10 mm
7411-0611*	Size 6	11 mm
7411-0612*	Size 6	12 mm
7411-0615*	Size 6	15 mm

* Special size only on request / optionally available

Implant material:

Highly crosslinked UHMWPE with vitamin E in accordance with ISO 5834-2

Implants

4-motion® insert VEX, highly crosslinked UHMWPE with vitamin E

Art.-no.	Size	Height
7411-0706	Size 7	6 mm
7411-0707*	Size 7	7 mm
7411-0708	Size 7	8 mm
7411-0709*	Size 7	9 mm
7401-0710	Size 7	10 mm
7411-0711*	Size 7	11 mm
7411-0712*	Size 7	12 mm
7411-0715*	Size 7	15 mm
7411-0806	Size 8	6 mm
7411-0807*	Size 8	7 mm
7411-0808	Size 8	8 mm
7411-0806*	Size 8	9 mm
7411-0810	Size 8	10 mm
7411-0811*	Size 8	11 mm
7411-0812*	Size 8	12 mm
7411-0815*	Size 8	15 mm
7411-0906	Size 9	6 mm
7411-0907*	Size 9	7 mm
7411-0908	Size 9	8 mm
7411-0909*	Size 9	9 mm
7411-0910	Size 9	10 mm
7411-0911*	Size 9	11 mm
7411-0912*	Size 9	12 mm
7411-0915*	Size 9	15 mm

Art.-no.	Size	Height
7411-1006	Size 10	6 mm
7411-1007*	Size 10	7 mm
7411-1008	Size 10	8 mm
7411-1009*	Size 10	9 mm
7401-1010	Size 10	10 mm
7411-1011*	Size 10	11 mm
7411-1012*	Size 10	12 mm
7411-1015*	Size 10	15 mm

4-motion® patella cemented

Art.-no.	Type	Dia.
4101-0826	DOM	26
4101-0829	DOM	29
4101-0832	DOM	32
4101-0835	DOM	35

* Special size only on request / optionally available

Implant material:

UHMWPE in accordance with ISO 5834-2

Implants

4-motion® PE insert (optionally available)

Art.-no.	Size	Height
7401-0106	Size 1	6 mm
7401-0108	Size 1	8 mm
7401-0110	Size 1	10 mm
7401-0112	Size 1	12 mm
7401-0115*	Size 1	15 mm
7401-0206	Size 2	6 mm
7401-0208	Size 2	8 mm
7401-0210	Size 2	10 mm
7401-0212	Size 2	12 mm
7401-0215*	Size 2	15 mm
7401-0306	Size 3	6 mm
7401-0308	Size 3	8 mm
7401-0310	Size 3	10 mm
7401-0312	Size 3	12 mm
7401-0315*	Size 3	15 mm
7401-0406	Size 4	6 mm
7401-0408	Size 4	8 mm
7401-0410	Size 4	10 mm
7401-0412	Size 4	12 mm
7401-0415*	Size 4	15 mm
7401-0506	Size 5	6 mm
7401-0508	Size 5	8 mm
7401-0510	Size 5	10 mm
7401-0512	Size 5	12 mm
7401-0515*	Size 5	15 mm

Art.-no.	Size	Height
7401-0606	Size 6	6mm
7401-0608	Size 6	8mm
7401-0610	Size 6	10mm
7401-0612	Size 6	12mm
7401-0615*	Size 6	15mm
7401-0706	Size 7	6mm
7401-0708	Size 7	8mm
7401-0710	Size 7	10mm
7401-0712	Size 7	12mm
7401-0715*	Size 7	15mm
7401-0806	Size 8	6mm
7401-0808	Size 8	8mm
7401-0810	Size 8	10mm
7401-0812	Size 8	12mm
7401-0815*	Size 8	15mm
7401-0906	Size 9	6mm
7401-0908	Size 9	8mm
7401-0910	Size 9	10mm
7401-0912	Size 9	12mm
7401-0915*	Size 9	15mm
7401-1006	Size 10	6mm
7401-1008	Size 10	8mm
7401-1010	Size 10	10mm
7401-1012	Size 10	12mm
7401-1015*	Size 10	15mm

* Special size only on request / optionally available

Implant material:

UHMWPE in accordance with ISO 5834-2

Instruments

4-motion® – General instrument set

No.	Name	Art. no.
1	Instrument Tray	9210000
2	Drillpin Adaptor 3.2	9210101
3	Drillpin 3.2, M3	9210201
4	Drillpin 3.2 with head	9210202
5	Pin Extractor	9210300
6	Cortical Palp 1.27 mm	9210401
7	Femoral Drill 6.4	9210501
8*	Femoral Trial 2s L	9210701
9	Femoral Trial 3s L	9210702
10	Femoral Trial 4s L	9210703
11	Femoral Trial 5s L	9210704
12	Femoral Trial 6s L	9210705
13	Femoral Trial 5w L	9210706
14	Femoral Trial 6w L	9210707
15	Femoral Trial 7w L	9210708
16*	Femoral Trial 8w L	9210709
17*	Femoral Trial 9w L	9210710
18*	Femoral Trial 2s R	9210801
19	Femoral Trial 3s R	9210802
20	Femoral Trial 4s R	9210803
21	Femoral Trial 5s R	9210804
22	Femoral Trial 6s R	9210805
23	Femoral Trial 5w R	9210806
24	Femoral Trial 6w R	9210807
25	Femoral Trial 7w R	9210808
26*	Femoral Trial 8w R	9210809
27*	Femoral Trial 9w R	9210810

No.	Name	Art. no.
28	Trial Insert A h=6	9211106
29	Trial Insert A h=8	9211108
30	Trial Insert A h=10	9211110
31	Trial Insert A h=12	9211112
32	Trial Insert B h=6	9211206
33	Trial Insert B h=8	9211208
34	Trial Insert B h=10	9211210
35	Trial Insert B h=12	9211212
36	Trial Insert C h=6	9211306
37	Trial Insert C h=8	9211308
38	Trial Insert C h=10	9211310
39	Trial Insert C h=12	9211312
40	Trochlea Rasp	9211410
41	Prime Handle	9211500
42	Attachment for Femoral / Tibial	9211501
43	Attachment for Insert Impactor	9211502
44	Attachment for Femoral Extractor	9211503
45	Alignment Handle	9211510
46	Alignment Rod	9211511
47	Alignment Rod Extension	9211512

NOTE:

The 4-motion® instrument set is split into two parts containing customised, patient-specific and system-specific, general instruments respectively.

For successful implantation of the system, the preoperative planning must be performed via the web-based platform operated by MEDivation and the patient-specific, single-use instruments ordered.

* Optional instrument

Instruments

4-motion® – General instrument set

No.	Name	Art. no.
63*	Ligament Balancer	9211600
64*	Torque Wrench for LB	9211601
65*	AO Adaptor for Torque Wrench	9211602
66*	Femoral Correction Block	9211700
68*	Adjustable Tibia Correction Bl. left	9211702
69*	Adjustable Tibia Correction Bl. right	9211703
70*	Angle Socket 0°	9211520

No.	Name	Art. no.
71*	Angle Socket 1°	9211521
72*	Angle Socket 2°	9211522
73*	Angle Socket 3°	9211523
74*	Distal Adjustable Femur Correction Bl.	9211704
75*	External Tibial Alignment Guide	9211001

4-motion® – Tibia PSA

No.	Name	Art. no.
48a	Drill and Keel Tower PSA	9211411
49a	Keel Punch PSA 1–3	9211412
50a	Keel Punch PSA 4–6	9211413
51a	Keel Punch PSA 7–10	9211414
52a	Tibial Drill PSA	9210503
53a*	Tibial Base Plate PSA 1	9210911
54a	Tibial Base Plate PSA 2	9210912
55a	Tibial Base Plate PSA 3	9210913
56a	Tibial Base Plate PSA 4	9210914
57a	Tibial Base Plate PSA 5	9210915
58a	Tibial Base Plate PSA 6	9210916
59a	Tibial Base Plate PSA 7	9210917
60a	Tibial Base Plate PSA 8	9210918
61a	Tibial Base Plate PSA 9	9210919
62a*	Tibial Base Plate PSA 10	9210920

4-motion® – Tibia 3° (optionally available)

No.	Name	Art. no.
48	Drill and Keel Tower	9211400
49	Keel Punch 1-3	9211401
50	Keel Punch 4-6	9211402
51	Keel Punch 7-10	9211403
52	Tibial Drill	9210502
53*	Tibial Base Plate 1	9210901
54	Tibial Base Plate 2	9210902
55	Tibial Base Plate 3	9210903
56	Tibial Base Plate 4	9210904
57	Tibial Base Plate 5	9210905
58	Tibial Base Plate 6	9210906
59	Tibial Base Plate 7	9210907
60	Tibial Base Plate 8	9210908
61	Tibial Base Plate 9	9210909
62*	Tibial Base Plate 10	9210910

* Optional instrument

Information:

Further information on the auxiliary instruments and their use is available from your customer representative, your distributor or directly from the manufacturer.



ARTIQO GmbH
Hans-Böckler-Straße 57
D-59348 Lüdinghausen / Germany

T +49 2591 / 89315-00
F +49 2591 / 89315-10
info@artiqo.de
www.artiqo.de

You can find the digital version
of the surgical technique at
[https://artiqo.de/download/
st-4-motion/](https://artiqo.de/download/st-4-motion/) or directly using this
QR code.

