
4-motion® Kniesystem

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4-motion®

artiqo.de

4-motion® Kniesystem

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DE	4-motion® Kniesystem	3-7
EN	4-motion® total knee replacement system	8-12

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1. General information

Device description

The 4-motion® total knee replacement system consists of three components that are designed to replace the femoral, tibial, and (optionally) the patellar articular surfaces of the knee joint.

Femoral and tibial base components:

The tibial base components and the CR femoral components containing the cruciate ligament are made of a Co-Cr-Mo alloy in accordance with ISO 5832-4, and the PSA tibial component is made of titanium alloy Ti6Al4V in accordance with ISO 5832-3. Models are available in different sizes and in left and right versions for cemented and non-cemented application. On request, the femoral component is also available in an allergy-inhibiting version with a titanium-niobium nitride (TiNbN) coating for allergy sufferers.

Tibial inlay component:

The "sliding surface" tibial inlays are part of the tibial components, and are made of highly crosslinked ultra-high-molecular-weight polyethylene "UHMWPE" with vitamin E ("VEX") in accordance with ISO 5834-2 and STM F2695-12. The tibial inlays are offered in various heights and sizes to fit the respective tibial base components.

Patellar component:

The single-component patellar component is made of ultra-high-molecular-weight polyethylene "UHMWPE" in accordance with ISO 5834-2 and is available in various diameters.

The 4-motion® knee endoprosthesis can be used with or without an artificial patella replacement, depending on the applied technique.

Contents of the packaging

- The device, contents of the packaging and materials used are defined by the device labels. This device must be employed in accordance with the respective, device-specific surgical technique.
- The lot and SN number(s) of the implants used must be documented in the patient's records. Corresponding labels are included with the packaging of the sterile implants for this purpose.
- The implant is supplied without bone cement.

Instructions regarding reuse

- The reuse of a previously inserted implant is prohibited.
- The implants are supplied sterile (sterilised using gamma radiation). They must not be re-sterilised.

Warnings/precautions

- The warnings on the packaging must be strictly observed.
- This device may only be used in its original condition. Modifications of any kind and mechanical manipulations are not permitted.

Information concerning materials used

The femoral component and the 3° tibial component are made of a cobalt-chrome alloy in accordance with ISO 5832-4, with the following chemical composition:

CoCrMo in accordance with ISO 5832-4	
Cr	26.5-30.0 wt%
Mo	4.5-7.0 wt%
Ni	Max. 1.0 wt%
Fe	Max. 1.0 wt%
C	Max. 0.35 wt%
Mn	Max. 1.0 wt%
Si	Max. 1.0 wt%
Co	Remainder

The PSA tibial component is made of a titanium-aluminium-vanadium alloy (Ti6Al4V) in accordance with ISO 5832-3 with the following chemical composition:

Ti6Al4V in accordance with ISO 5832-3	
Al	5.50-6.75 wt%
V	3.50-4.50 wt%
Fe	Max. 0.30 wt%
O ₂	Max. 0.20 wt%
C	Max. 0.08 wt%
N ₂	Max. 0.05 wt%
H ₂	Max. 0.015 wt%
Ti	Remainder

The femoral and tibial components for non-cemented anchoring have a special titanium coating (titanium plasma spray, TPS, grade 3 pure titanium in accordance with ISO 5832-2) on the posterior surfaces. The titanium plasma spray coating in accordance with ISO 5832-2 corresponds to the following chemical composition:

TPS grade 3 in accordance with ISO 5832-2	
N	0.05 wt%
C	0.08 wt%
H	0.0125 wt%*
Fe	0.30 wt%
O	0.35 wt%
Ti	Remainder

*Except for billets, for which the maximum hydrogen content shall be 0.010% (mass fraction) and for flat products for which the maximum hydrogen content shall be 0.015% (mass fraction).

For patients who are known to have a metal intolerance, there is also a special "sensitive version" of the femoral component with a coating made of titanium-niobium nitride (TiNbN).

The patellar component is made of ultra-high-molecular-weight polyethylene (UHMWPE, GUR 1020 type 1), a long-chain-hydrocarbon polymer. The used material complies with the qualities of ISO 5834-1, ISO 5834-2 and ASTM F648. In accordance with standard ISO 5834-1, the permissible impurities per kg in the powder material are below the following limit values:

Maximum permissible impurities in UHMWPE in accordance with ISO 5834-1	
Ash	125 mg/kg
Ti	40 mg/kg
Ca	5 mg/kg
Cl	30 mg/kg
Al	20 mg/kg

The highly crosslinked VEX inlay with vitamin E is made of the same base material and contains 0.1 percent by weight of α -tocopherol (vitamin E). The material for the VEX inlay meets the requirements of ASTM F2695-12.

All materials used are clinically proven and correspond to the state of the art for knee joint implants.

2. Storage and handling

Storing the sterile implants

Implants should always be stored in their unopened original packaging. They should be stored clean and dry, protected from direct sunlight and at an ambient temperature between 10°C and 30°C.

The materials have been chosen in such a way that temporary temperature fluctuations, as can occur during transport, have no negative effects on the performance of the implants and their packaging.

Handling the sterile implants

The cardboard box is sealed with an ARTIQO Security Seal. If the seal label displays abnormalities (e.g. defects or tears which no longer guarantee that the packaging is sealed or the white ARTIQO logo on the red background is not clearly visible and/or the Security Seal has been broken ("Void"/"Open" is legible in such cases)), it is no longer possible to guarantee the integrity of the protective packaging.

Prior to insertion of the sterile implant, the packaging must be inspected for damage. If the primary packaging is damaged, the sterility is no longer guaranteed and the implant must not be used. Should it be evident that air has penetrated the sterile barrier system (not airtight) or the sterilisation indicator not be red, the sterility is also no longer guaranteed and the device must not be used.

The implants are vacuum-packed in a triple peel pouch (sterile barrier system). The packaging system consists of two inner sterile pouches made of composite film that can be peeled back to open the pouches. A third outer pouch with the same specifications serves as protective packaging for the device. The middle (second) pouch features a label including a symbol for the structure of the sterile barrier system, as well as a red indicator dot.

When unpacking the implant, a check must be performed to verify that it matches the designation on the packaging (art. no. and size).

The corresponding aseptic regulations for surgical personnel must be observed when removing the implant from the packaging; the outer and, if appropriate, central PE pouch are opened by the non-aseptic surgical assistant; the inner pouch is opened by the aseptic surgical assistant. The implant must not come into contact with objects which could damage its surface. Each implant must be visually inspected for damage prior to insertion.

		Inlay									
		Femur first									
		1	2	3	4	5	6	7	8	9	10
Patella	Ø 26										
	Ø 28										
	Ø 32										
	Ø 35										

Tibia / Tibia PSA

3. Instructions for use – special application instructions

Patient briefing

Each patient must be comprehensively briefed on the procedure and possible risks prior to the surgery.

Patients receiving a knee replacement should be made aware that the useful service life of the implant depends on their weight and level of activity. Excessive physical activity or trauma can result in loosening, disproportionate wear and/or breakage of the implant. The patient must be informed of the implant's capabilities and possible consequences for their lifestyle. The implant may not last for the rest of the patient's natural life or a specific period. Implants are not as resilient and durable as natural, healthy tissue and bones. As such, it is possible that all components will need to be replaced after a certain period of time.

Following the implantation, the patient is to be provided with a completed Implant Card, which allows traceability of all implanted components.

Postoperative care must be performed in accordance with the applicable guidelines and hospital regulations in coordination with the patient.

Preoperative planning

The preoperative planning provides essential information for the selection of the appropriate components and possible component combinations. The planning of the 4-motion® total knee replacement system can be performed either conventionally using X-ray templates (digital or analogue). In this case, the surgical procedure is performed using reusable, conventional cutting blocks. Alternatively, planning can be carried out using a three-dimensional model of the patient's pathologically altered knee via the web-based planning platform provided by Medivation. This planning is based on the patient's CT data. Based on this planning, patient-specific single-use cutting blocks are manufactured, which are adapted to the individual surface contour and anatomy of the respective knee joint.

Surgical technique

The implantation must be performed in accordance with the valid version of the corresponding surgical technique. The latest version is available to download on the ARTIQO GmbH website (www.artiqo.de/en/downloads). Familiarity with the recommended surgical technique for this system and its careful performance is essential for optimal results.

Combination restrictions

The 4-motion® total knee replacement system is designed as a complete system and does not allow for any combination or exchange with components of other systems or manufacturers.

The "PE sliding surface" tibial inlays are designed to be the same size as the tibial base component.

The femoral component should be no more than two sizes different from the tibial component. The patellar components can be combined with all femoral component sizes.

The cement manufacturer's instructions for use must be observed when using bone cement.

The instrument set supplied by ARTIQO is coordinated to the implants; other instruments must not be used. (With the exception of standard surgical instruments.)

4. Application instructions

General risks of a surgical procedure are not described in these instructions for use. The surgeon must be proficient in the recognised surgical techniques, both in theory and in practice, and is responsible for the proper execution of the surgical procedure.

Damaged or surgically removed implant components must not be reused. Implants that have been used once must not be reused. There may be defects that are invisible to the naked eye that could reduce the lifetime of the implant.

The preoperative assessment of the patient's anatomical suitability is carried out based on clinical and radiological examinations and/or CT scans. Only patients who meet the criteria listed in the section on indications for use should be selected for the surgery.

Additional implants or sizes should be available as a fallback option in case the intended implant cannot be used.

If a metal allergy or intolerance is suspected, the patient should be given the implant with special coating.

It is important that the surgeon and surgical team are familiar with the system-specific surgical technique, the implant range and the instrument set, and that all of the necessary system documentation is available on site in full.

Implant instruments, including any planned or ordered patient-specific single-use cutting-blocks, must be available and fully functional prior to the procedure.

The cement manufacturer's instructions for use must be observed when using bone cement. The implant site must be flushed prior to insertion of the cement (in cemented anchoring). It must be ensured that all loose particles (e.g. bone fragments) are removed from the prepared implant site. Once cementing is complete, all excess or loose cement particles must be removed from the wound site.

The artificial joint replacement may become loose due, for example, to overuse, wear, poor bone quality or infection.

The revision of a knee prosthesis is a complicated joint replacement procedure. It should always be expected that revision surgery will be required. In this case, a ready-for-use revision system should be available on site.

Indications

- Patients with advanced degeneration of the knee joint resulting from degenerative and post-traumatic arthritis (PTA) or rheumatoid arthritis (RA)
- Avascular necrosis of the femoral condyle
- Moderate varus, valgus or flexion deformity

Contraindications

- Acute or chronic infections, local or systemic
- Severe muscular, neural or vascular diseases putting the affected extremities at risk
- Absent bone substance, bone tumours or inadequate bone quality posing a risk to the stable fit of the prosthesis
- Severe osteoporosis or osteomalacia
- Hypersensitivity to the material employed
- Loss of the ligaments
- Any concomitant disease capable of putting the function of the implant at risk
- Revision with pronounced bony defects

Precautions/restriction

- Major malpositions in combination with instabilities
- BMI > 30
- Loss of posterior cruciate ligament
- High level of physical activity (e.g. competitive sports, strenuous physical work)
- Diabetes

Possible side effects and complications

The adverse effects listed below are some of the most typical and common consequences of a knee arthroplasty.

- Loosening, wear and breakage of the implant components
- Early- or late-stage infections rendering removal of the implant or more extensive surgical procedures necessary.
- Pain

- Dislocation, partial dislocation
- Restricted range of motion
- Restricted functionality
- Instability
- Fractures of the tibia or femur
- Joint effusion
- Arthrofibrosis
- Necrosis of the patella
- Soft tissue impingement (e.g. "patellar clunk")
- Delayed healing or wound dehiscence
- Haematoma
- Nerve or vessel damage
- Local intolerances and tissue reactions
- Toxic systemic side effects
- Allergies to material components
- Venous thrombosis and pulmonary embolism
- Cardiovascular and pulmonary disorders (e.g. fat embolism syndrome)
- Complications in connection with blood transfusions
- Change in position (migration) of one or more prosthesis components
- Differences in leg length
- Excessive wear of the polyethylene components due to intraoperative damage of the femoral component, foreign bodies such as loose cement and/or bone fragments and/or the patient having an excessive level of activity or excessive body weight.
- Periprosthetic osteolysis
- Tendinitis in the operated extremity
- Periarticular calcification or ossification, with or without impact on the mobility of the joint.

All these complications must be diagnosed and treated as early as possible and may require revision surgery.

Risk factors for the side effects and complications listed above are determined by patient-specific factors such as level of activity, weight, neuromuscular status and concomitant diseases, traumatic events as well as surgical factors such as size selection, positioning of the components and/or alignment or cementing technique.

The following factors may be detrimental to the success of the procedure:

- Osteoporosis
- Deformations
- Local bone tumours
- Systemic diseases and metabolic disorders
- Reference to infections and falls in the medical history
- Drug addiction or abuse
- Obesity
- Major physical activities and activities associated with strong vibrations, where the prosthesis is subjected to impacts and/or excessive loading (e.g. strenuous physical work, long-distance running, etc.).

Avoidance/reduction of side effects

Important: If the implantation of this knee prosthesis is deemed the best solution for the patient and some of the circumstances described above apply in the case of the patient, it is particularly important to brief the patient about the expected consequences of these circumstances for the success of the surgery.

Additionally, it is recommended that the patient be informed of activities which can be utilised to reduce the consequences of these complicating circumstances. All information provided to the patient must be documented in writing by the surgeon performing the procedure.

The surgeon must be familiar with the anchoring philosophy of the implants and the surgical technique.

Proper selection of the implant size and position reduces the likelihood of post-operative migration.

The use of antibiotic-loaded bone cement can reduce the risk of prosthetic joint infection (PJI) in primary total joint arthroplasty (TJA).^[1] When high-viscosity bone cement is used, a correlation with an increased risk of aseptic loosening has been observed in the literature with different systems.^[2] When polyethylene inlays with a height of ≥ 13 mm are used, a correlation with an increased risk of aseptic loosening has been observed with CR knee replacement systems.^[3]

Postoperative treatment

Recognised procedures should be employed for the postoperative care and treatment. The documentation of the postoperative treatment should be performed in accordance with the instructions and regulations in force at the hospital. For example:

- Surgical report
- Postoperative x-rays
- Patient Implant Card
- Regular check-ups/follow-up examinations

Interactions

The 4-motion® knee prosthesis contains paramagnetic materials and thus only reacts weakly to strong magnetic fields.^[4] The higher the field strengths, the more pronounced the interactions. Reactions and effects due to heating, migration and image artefacts can be expected during MRI scans of patients with implanted endoprostheses.

The patients must be made aware that extremely strong electromagnetic fields are to be avoided and the technical specifications in the MRI safety information must be taken into consideration for scans. Furthermore, an artificial knee joint may set off metal detectors (e.g. at airport security).

MRI safety information

Non-clinical testing in accordance with the US FDA Directive and the ASTM standards F2052, F2182, F2213, and F2119 has confirmed that the TKA from ARTIQO can be classed as “MR conditional”. Consequently, the products have been labelled in accordance with ASTM F2503.

Patients with an ARTIQO 4-motion® knee prosthesis according to the list of implant components tested for the 4-motion® total knee replacement system can be scanned safely under the tested conditions specified in the MRI safety information for the 4-motion® total knee replacement system.

Failure to comply with these conditions may result in patient injury.

Further information on MRI safety can be found at <https://artiqo.de/en/downloads>. The MRI safety information for the 4-motion® knee system can be requested there.

Disposal of non-compliant implants and handling of explants

Non-compliant, defective implants (lapsed use-by date, defective packaging, etc.) should be returned to the manufacturer or professionally rendered unusable and disposed of by the owner.

Implant revisions should be performed in accordance with the state of the art of science and technology as well as the specific implant based on the surgical technique.

Explants should be treated as infectious waste and disposed of or reprocessed in accordance with the hospital regulations.

Please note that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority without delay.

^[1] Sebastian et al. 2020: Antibiotic containing bone cement in prevention of hip and knee prosthetic joint infections: A systematic review and meta-analysis. J Orthop Translat. 2020 May 8;23:53-60. doi: 10.1016/j.jot.2020.04.005.

^[2] Buller et al. 2020: Primary Total Knee Arthroplasty Performed Using High-Viscosity Cement is Associated With Higher Odds of Revision for Aseptic Loosening. J Arthroplasty. 2020 Jun;35(6S):S182-S189. doi: 10.1016/j.arth.2019.08.023.

^[3] Rajamäki et al. 2022: Thicker polyethylene inserts (≥ 13 mm) increase the risk for early failure after primary cruciate-retaining total knee arthroplasty (TKA): a single-centre study of 7643 TKAs. Knee Surg Sports Traumatol Arthrosc. 2023 Mar;31(3):1018-1025. doi: 10.1007/s00167-022-07189-8.

^[4] Frank G Shellock: reference manual for magnetic resonance safety, implants and devices. (Yearly updated). Detailinformationen: <http://www.mrisafety.com/>

5. Handling, care and reprocessing of instruments

Section 5 applies only for instruments – implants must never be reprocessed.

The following recommendations are intended for information purposes only. The manufacturer cannot be held liable for any instrument cleaning and sterilisation processes performed at the purchaser's facilities.

The instrument set is supplied “non-sterile” and must be cleaned, disinfected and sterilised using validated procedures prior to being used.

Brief summary for reprocessing instruments:

Our instructions are based on validation processes performed as part of the preparation of the medical device for reuse. The reprocessor is responsible for ensuring that the reprocessing actually performed with the equipment and personnel available in the reprocessing facility achieves the required results.

Precleaning:

Prior to every automated cleaning, precleaning should be performed as quickly as possible after use and preferably using an ultrasound bath. Some washer disinfectors are not capable of reprocessing instruments without residues, especially after surgical residues have been left to dry for extended periods of time.

Automated cleaning and disinfection

The automated cleaning and disinfection should be performed using a process validated in accordance with the applicable rules and standards in a washer disinfectant in accordance with EN ISO 15883-1 and -2, preferably using (mild) alkaline cleaning agents.

Sterilisation

The steam sterilisation is performed in a steam steriliser using a process validated in accordance with EN ISO 17665 employing the fractionated vacuum process as per EN 285. The prevalidation is performed using the following programme specifications:

Cycle type	Fractionated vacuum
Temperature	134 °C *
Pressure	3050 mbar abs.
Sterilisation time	5 minutes
Drying time	25 minutes

*Sterilisation temperature of 134°C (273°F) plus tolerance as per ISO 17665. Please observe the respective nationally applicable regulations.

Other sterilisation methods and cycles may also be used. However, they must be validated in accordance with the applicable rules and standards before being adopted into routine use.

Inspection, repairs and testing:

Dry wet or damp devices.

Clean and disinfect any visibly contaminated devices again.

After every cleaning, disinfection and drying, check that the instruments are dry, clean and free from damage.

Some instruments feature colour marks made of ceramic, which degrade over time with repeated reprocessing. Check the colour marks for visible or tangible changes (cracks, detachments, loosening, protrusions, etc.)

If proper functioning is no longer guaranteed, wear is evident on the instrument or a change to the colour mark is identified (cracked, loose, slightly protruding colour mark), the items in question must be repaired or replaced.

Assemble any instruments that can be dismantled. Only instruments which have been checked and are functioning correctly and ready for use should be stored. Check the compatibility with the corresponding devices (couplings).

Packaging and transport


Please use our transport packaging for returns and include our cleaning and sterilisation form.

For more detailed information on the preparation, precleaning, disinfection and sterilisation/resterilisation, inspection and care, packaging and storage of our instruments, please consult our reprocessing instructions (NA30-04-05).

6. Further information:

The current Summary of Safety and Clinical Performance (SSCP) of the 4-motion® total knee replacement system can be requested via the website www.artiqo.de/en/downloads. For further information on the use of this device, please contact your local sales office.

7. Explanation of the symbols used by ARTIQO GmbH

	Do not re-use
	Do not re-sterilize
	Non-sterile
	Sterilized using irradiation
	Double sterile barrier system with outer protective packaging
	Indicator for radiation sterilization (changes from yellow to red/violet following gamma irradiation)
	Medical device
	Batch code
	Serial number
	Unique Device Identifier
	Manufacturer
	Do not use if package is damaged
	Keep away from sunlight
	Keep dry
	CE marking (if applicable, with 4-digit number of the notified body)
	Temperature range for storage temperature
	Caution
	Consult instructions for use
	Date of manufacture
	Use-by date
	Catalogue number
	MR conditional
	Quantity

ARTIQO GmbH
Hans-Böckler-Straße 57
59348 Lüdinghausen
Germany

T +49 2591 / 89315-00
F +49 2591 / 89315-10
E-Mail: info@artiqo.de
www.artiqo.de

eIFU

Download:

www.artiqo.de/downloads/ifu