ELEC[®] and ELEC[®]plus hip implants

Purity + Strength = Permanence

Instruction for use

HiPer Medical AG

Ziegeleistrasse 7
16727 Oberkrämer OT Eichstädt
Germany
hiper@hiper.de
www.hiper-medical.de
Tel.: +49 (0) 3304 3977-0
Fax.: +49 (0) 3304 3977-90

DE	Gebrauchsanweisung ELEC® und ELEC®plus Implantate
BE	Импланти ELEC® и ELEC®plus Инструкция за използване
BG	Импланти ELEC® и ELEC®plus Инструкция па выкарыстанні
CS	Návod k použití Implantáty ELEC® a ELEC®plus
DA	Brugsanvisning ELEC®- og ELEC®plus-implantater
EL	Οδηγίες χρήσης Εμφυτεύματα ELEC® και ELEC®plus
EN	Instructions for use ELEC® und ELEC®plus Implants
ES	Instrucciones de uso Implantes ELEC® y ELEC®plus
ET	Kasutusjuhendit Implantaadid ELEC® ja ELEC®plus
FI	Käyttöohjetta ELEC®- ja ELEC®plus-implantit
FR	Mode d'emploi Implants ELEC® et ELEC®plus
HR	Upute za korištenje Implantati ELEC® i ELEC®plus
HU	Használati utasítás ELEC®- és ELEC®plus-implantátumok
IT	Manuale d'istruzioni Impianti ELEC® e ELEC®plus
LT	Naudojimo instrukcijos ELEC® ir „ELEC®plus“ implantai
LV	Lietošanas instrukcija ELEC® un ELEC®plus implanti
NL	Gebruiksaanwijzing ELEC®- en ELEC®plus-implantaten
NO	Bruksanvisning ELEC®- og ELEC®plus-implantater
PL	Instrukcja użycia Implanty ELEC® i ELEC®plus
PT	Instruções de utilização implantes ELEC® e ELEC®plus
RO	Instructiuni de utilizare Implanturi ELEC® și ELEC®plus
RU	Инструкция по применению имплантатов ELEC® и ELEC®plus
SK	Návod na použitie Implantáty ELEC® a ELEC®plus
SI	Navodila za uporabo Vsadki ELEC® in ELEC®plus
SV	Bruksanvisning ELEC® och ELEC®plus implantat
TR	Kullanım Talimatı ELEC® ve ELEC®plus implantları
UK	Імпланти ELEC® та ELEC®plus Інструкція із застосування
VI	Hướng dẫn sử dụng Implants ELEC® và ELEC®plus
ZH	ELEC® 和 ELEC®plus 植入物

DE – Deutsch

ELEC®- und ELEC®plus-Implantate

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

Dobbelt sterilt barriersystem

**ELEC® and ELEC®plus implants**

Before using these devices, the user is obliged to carefully study the following recommendations and instructions, together with the device-specific instructions. The distributor of these devices does not accept any liability for direct or subsequent damage resulting from improper use or handling, in particular due to non-compliance with the following instructions for use, or due to inappropriate maintenance or servicing. These implants may only be used by doctors with adequate experience and expertise in hip replacement surgery. Familiarity with the surgical technique recommended for this system and its careful application are vital to achieve the best possible outcome.

An ELEC® or ELEC®plus femoral head may only be paired with an ELEC®plus acetabular insert or a PE cup approved by the manufacturer of the other components or an approved PE acetabular insert.

An ELEC®plus acetabular insert may only be paired with an ELEC® or an ELEC®plus femoral head or ELEC®plus revision femoral head. Pairing with femoral heads from another manufacturer is prohibited.

The instructions for use from the manufacturer of the other components (endoprosthesis manufacturer) must be observed.

1 Product description and implant raw materials

ELEC® and ELEC®plus femoral heads should only be used for fixing to stems approved by the endoprosthesis manufacturer.

ELEC®plus acetabular inserts are only intended for fixing in acetabular housings approved by the endoprosthesis manufacturer.

ELEC® and ELEC®plus are high-purity implant ceramics (ISO 6474-1 and ISO 6474-2). These raw materials are biologically inert, biocompatible, biostable, mechanically stable, corrosion-resistant and they have excellent resistance to fatigue, high shock resistance and tensile strength, excellent fracture toughness, an extremely high degree of hardness and they do not release ions into the body. The device, package content and used raw materials are defined on the device labels.

2 Handling**2.1 General instructions**

Implants must always be stored in their complete, unopened protective packaging. The packaging of the implants should not be exposed to direct sunlight. Before inserting the implant, the packaging should be checked for damage, as this can impair the sterility of the implant. When unpacking the implant, it should be checked that it corresponds to the description on the packaging (serial number/size). When removing the implant from the packaging, appropriate hygiene standards must be observed (see point 3). It must be ensured that all the surfaces of the implant are protected from damage, as this could be a determining factor in the possible failure of the implant. Therefore, the implant should not come into contact with objects that could damage its surface. Each implant should undergo a visual inspection for damaged areas prior to use. Altering or bending an implant may not only shorten its lifetime, but also result in the immediate or subsequent failure of the prosthesis when it is loaded. The implant must therefore not be altered mechanically or in any other manner. Implants from damaged packaging, as well as non-sterile, contaminated, damaged, improperly handled implants, implants that have been altered without authorisation or implants that have fallen on the floor must not be used.

Caution ELEC® and ELEC®plus implants are intended for single use only. It is prohibited to reuse a prosthesis component that had previously been implanted in the body. Reuse may lead to breakage due to scratches or impact marks.

2.2 Instructions for use

Only the system instruments recommended by the manufacturer of the other components should be used for the implantation. The enclosed instructions for use should be noted

before the instruments are used. ELEC® and ELEC®plus implants may only be used in the combinations recommended by the manufacturer of the other components.

Only ELEC® and ELEC®plus implants with taper dimensions that match those of the combination implant may be used. The taper size is indicated on the label and, if possible, also on the implant. Only brand new implants should be used.

ELEC® and ELEC®plus femoral heads

The tapered plug-in connection enables reliable clamping of the modular femoral head. It brings about an even transfer of forces to the stem. The ELEC® and ELEC®plus femoral heads must fit perfectly with the taper of the prosthesis stem.

ELEC®plus acetabular inserts

The clamping conical connection enables a reliable mounting of the modular acetabular insert in the implanted acetabulum. It brings about an even transfer of forces into the cup. The ELEC®plus acetabular inserts must fit perfectly with the implanted cup. With press-fit fixation in particular, the cup may be slightly deformed following implantation in the bone bed, which can hinder the correct fixation of an ELEC®plus acetabular insert.

2.3 Surgical technique

2.3.1 ELEC® and ELEC®plus femoral heads

Attaching the femoral head to the prosthesis stem

In order to guarantee reliable function, it is essential that the utmost care be taken. The plastic cap protecting the stem taper from damage should only be removed just before attaching the test head. The following tasks should be completed before attaching the ELEC® or ELEC®plus femoral head:

1. Rinse the stem taper thoroughly with water to remove any tissue, bone fragments or cement residue
2. Dry the stem taper with a clean gauze pad and check the position and function of the joint after attaching the test head.
3. Remove the test head and repeat the rinsing and drying steps
The following tasks should be completed when attaching the ELEC® or ELEC®plus femoral head:
4. Place the ELEC® or ELEC®plus femoral head onto the stem taper with a little axial pressure and simultaneous rotation until it is firmly in place
5. Now place the plastic impactor on the pole of the head and finally clamp the femoral head onto the stem with a light hammer blow in the axial direction. Hitting with the hammer causes plastic deformation of the surface structure of the stem taper and optimal pressure distribution is achieved. This leads to non-twisting fixation.

Caution Never hit a ceramic femoral head with a metal hammer.

Damaged stem taper

An ELEC® or ELEC®plus femoral head that has been hammered once onto a stem taper and removed again must not be reattached. A metal head or ELEC®plus revision femoral head must be used instead. The same applies if the prosthesis stem is left in situ during a revision, but the femoral head is replaced.

If an ELEC®plus acetabular insert was intended as the articulating component in any of the above cases, it may have to be exchanged for a PE acetabular insert. The specifications of the endoprosthesis manufacturer must be observed.

Caution It is only possible to leave the stem in situ if an inspection finds that the stem taper is flawless.

2.3.2 ELEC®plus acetabular inserts

Securing the acetabular insert in the implanted acetabulum.

In order to guarantee reliable function, it is essential that the utmost care be taken. The following tasks should be completed before the ELEC®plus acetabular insert is inserted:

1. Ensure that the inclination (40° – 45°) and anteversion (10° – 20°) of the acetabulum are appropriate for the indication

2. Rinse the clamping cone of the implanted acetabulum thoroughly with water to remove any tissue, bone fragments or cement residue
3. Dry the clamping cone of the implanted acetabulum and check the position and function of the joint by inserting a trial implant. If any screws are used, they must be adequately deeply recessed.
4. Remove the trial implant and repeat the rinsing and drying steps

The following tasks should be completed if an ELEC[®]plus acetabular insert is used:

5. Grip part of the edge of the ELEC[®]plus acetabular insert between your thumb and forefinger and carefully slide it into the clamping cone of the acetabulum. Do not press down.
6. Ensure that the flush fit of the ELEC[®]plus acetabular insert is correct by feeling the edge of the cup with your finger. You should be able to feel that the flat surface of the cup and the flat surface of the ELEC[®]plus acetabular insert are evenly closed and, if required, you can correct the fit through careful, targeted pressure.
7. Only now, firmly press the ELEC[®]plus acetabular insert with your thumb in the centre of the cap and secure it in the clamping cone, and check again that there is a flush fit and - if provided - use the instrument approved by the manufacturer of the endoprosthesis system.

2.3.3 Reoperation due to a broken femoral head

In the event of a broken ceramic femoral head, the below procedure should be followed:

Remove all ceramic particles. A PE cup should be removed, even if it is firmly fitted because pushed-in ceramic particles may damage the new femoral head, which would again lead to increased PE wear. An ELEC[®]plus acetabular insert must also be removed if a metal femoral head is to be used. Under no circumstances should another ceramic femoral head be attached to the stem taper. A metal femoral head or a ceramic revision femoral head approved by the endoprosthesis manufacturer must be used. The prosthesis stem may only be left in situ if a thorough inspection finds that the stem taper is flawless.

2.3.4 Reoperation due to a broken ceramic acetabular insert

In the event of a broken ceramic acetabular insert, the below procedure should be followed:

Remove all ceramic particles. If it is not possible to remove all the ceramic particles with the cup still in situ or if the cup has been damaged due to the breakage of the ceramic acetabular insert, the cup must be removed and replaced by a new one after all the ceramic particles have been removed.

3 Packaging and sterility

The implants are delivered sterile (gamma sterilisation, at least 25 kGy) in double tear-open packaging. The specified expiry date is based on the packaging being undamaged and unopened, and stored according to the specified storage conditions. Sterilisation by gamma radiation may alter the colour of ceramic devices. This does not affect the stability or any other property of the implant.

Caution The implants must not be resterilised. It is not permitted to reprocess unimplanted components if their packaging has been opened because individual validated processes have to be run through again.

The implants are delivered in a double sterile barrier system. This enables an aseptic

presentation technique in which the non-sterile surgical personnel open the first bag and hand it over to the sterile surgical personnel to take out the second bag. The sterile surgical personnel can remove the second bag and place it on the sterile work station.

4 Preoperative planning

Preoperative planning using X-rays, CT data and other information is indispensable and provides important information about the suitable prosthesis type, placement and possible component combinations. The surgery should only be carried out if the patient's material compatibility has been positively determined. When planning the surgery, the X-ray templates of the manufacturer of the other components should be used. Test prostheses to

verify that the position is correct (as applicable) and additional implants should be available, in case other sizes are required or the intended implant cannot be used. Accepted surgical practices should be applied during postoperative care.

5 Indications

- Total hip replacement, primary surgery in combination with prosthesis stems and cups
- Revision (= endoprosthesis removal) using new prosthesis stems and cups
- Hemiarthroplasty (partial hip replacement, only applicable to femoral heads)

6 Contraindications

- Revision with retention of the stem and/or cup in situ (risk of component breakage due to damaged fastening cones)
- Only acetabular inserts: Inclination outside 40° – 45°, anteversion outside 10° - 20° (risk of luxation, subluxation, impingement)
- Only acetabular inserts: Retroversion

7 Possible adverse effects and complications

The negative effects listed below are among the most common consequences of hip replacement surgery:

- Dislocation, subluxation, restricted range of motion, unwanted change in leg length
- Infection, hematoma, delayed wound healing
- Venous thrombosis and pulmonary embolism and other cardiovascular disorders
- Displacement or loosening of the implant
- Fracture of the implant, bone or cement
- Side effects due to joint surface friction and wear (e.g. osteolysis, inflammatory tissue reactions)

8 Patient information, documentation

The serial numbers of the used implants should be documented in the patient records. Appropriate labels are enclosed with the packaging of the sterile implants. The advantages and risks of the procedure should be explained to the patient. If the implant is considered the best solution for the patient, even though the above-described contraindications partly apply to the patient, the patient must be informed about the expected effects of these circumstances and the expected risks. Patients who receive a hip replacement must be informed that the lifetime of the implant depends on their weight and level of activity. It is also recommended to inform the patient about activities that can help to reduce the impact of these aggravating circumstances. All information given to the patient must be documented in writing by the operating doctor. After the surgery, the patient should be given the enclosed implant card, which contains all the necessary information about the implant. Adhesive labels are provided for the documentation of the used implant. It should also be recorded on the implant card that future MRI examinations could be considered on a case-by-case basis due to possible interactions with the implant. The patient should also be informed about this.

9 Statement on MRI examinations

The ELEC® and ELEC®plus implants can be considered non-critical with regard to MRI examinations. The interactions of the ELEC® and ELEC®plus materials with magnetic fields are most likely insignificant due to their diamagnetic properties. However, the treating doctor is obligated to monitor the patient and the implant during MRI examinations to detect any unexpected adverse effects.

10 Summary of Safety and Clinical Performance

The Summary of Safety and Clinical Performance in accordance with Article 32 of the Medical Device Regulation 2017/745 is in the EUDAMED database.

SRN (single registration number) HiPer Medical AG: DE-MF-000008894

11 Explanation of the label symbols

Caution, consult accompanying documents

Date of manufacture (month/year)


Sterilized using irradiation

Do not resterilise

Serial number


Permissible storage temperature range

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