CEMSTOP®

REF	T770008 / T770010 / T770012 / T770014 / T770016 / T770018
-----	---

Basic UDI-DI / IUD-ID de base / UDI-DI básico / Basis-UDI-DI : 376017704B02CU

en	Bio-absorbable cement restrictor
	<i>Instructions for use</i>
fr	Obturbateur diaphysaire biorésorbable
	<i>Notice d'utilisation</i>
es	Restrictor de cemento bioabsorbible
	<i>Manual de instrucciones</i>
de	Bioresorbierbarer Zement-Restriktor
	<i>Gebrauchsanweisung</i>

Manufactured in France by / Fabriqué en France par / Fabricado en Francia por / In Frankreich hergestellt von :
TEKNIMED S.A.S
8, rue du Corps Franc-Pommies
65500 VIC EN BIGORRE - FRANCE
 Tél. (33)5 34 25 10 60
 www.teknimed.com

First CE marking date applied: 1998

GENERAL INFORMATION
 Before using medical devices, the operator must be thoroughly familiarised with the safety rules mentioned in the current instructions for use, as well as the information relating to each product (description, leaflets, etc.). The corresponding information is available from supplier. The operator should also be informed of the residual risk(s) associated with the product used.

Summary of Safety and Clinical Performance is available on EUDAMED database (https://ec.europa.eu/tools/eudamed) or https://myimplantcard.com/sscp/

INTENDED USERS
 Medical devices should only be used in an adapted environment and by healthcare qualified operators having a sound knowledge and full mastery in preparation techniques specific to such products. Preparation techniques may be acquired from qualified distributors. The operator is responsible for any complications or harmful consequences which might result from an erroneous indication or operative technique, an improper use of the equipment and/or a failure to comply with the safety rules provided in the instructions for use. Neither manufacturer nor qualified distributor can be held responsible for these complications.

INTENDED USE
CEMSTOP® cement restrictor is a diaphyseal plug intended to occlude the medullary cavity before the introduction of acrylic cement during hip arthroplasty surgeries. **CEMSTOP®** is available in the following six sizes: 8-10 - 12 - 14 - 16 - 18 mm. Each size corresponds to the diameter of the medullary canal to occlude.

LIFETIME
CEMSTOP® is an implantable device, intended to be totally dissolved and metabolised by the organism within a time lag of 1 to 15 days. This period may vary depending on the physiological condition of the patient. The degradation products of **CEMSTOP®** are gelatin, glycerin, water and methylparahydroxybenzoate.

INDICATIONS
CEMSTOP® cement restrictor is indicated in cemented hip arthroplasties.

TARGET POPULATION
CEMSTOP® is to be used in patients for whom a cemented hip arthroplasty surgery is indicated.

PERFORMANCES AND EXPECTED CLINICAL BENEFITS
CEMSTOP® prevents the cement from flowing down the diaphysis and therefore facilitates cement pressurization. It resists the temporary increase in temperature generated by the polymerization of acrylic cement.

CONTRAINDICATIONS
 - Procedures other than those stated in the INDICATIONS section.
 - Patients susceptible to allergic reactions to the **CEMSTOP®** absorbable components and their products of metabolism.
 - As of today, no information concerning possible contraindications of **CEMSTOP®** on children, pregnant or breastfeeding women are published. Nevertheless, this population remains at risk for all complications related to a surgical procedure.

SIDE EFFECTS AND ADVERSE EVENTS
 No undesirable side-effects have ever been recorded. Interaction with other agents: none known to date.

Any serious incident occurring in relation to the device **CEMSTOP®**

must be reported without delay to the manufacturer and the competent local authority where the user and/or the patient is established.

COMPOSITION	
Gelatin (porcine-based)	40% (w/w)
Glycerin	35% (w/w)
Water	25% (w/w)
Methylparahydroxybenzoate	2,5% (w/w)

STERILISATION
CEMSTOP® is delivered sterile under double packaging, ready to use in the operating room. The product is sterilized by gamma rays at 25kGy. Sterility is guaranteed until the expiry date if the packaging has not been opened or damaged.

RECOMMENDED ACCESSORIES
 - Complete set instrumentation (T067750) sold separately.

INSTRUCTIONS FOR USE
 The diameter of the cement restrictor can be predetermined by the surgeon by setting up the correspondence between the supplied templates of the prosthesis and diameters of the medullary canal.

- Select the sizing trial device corresponding to the medullar canal diameter predetermined with the templates or the rasps of the prosthesis ancillary, screw it on the inserter and check that it is correctly secured.
- Lock the set screw.
- Introduce the sizing trial device into the medullar canal to check for a proper fit. Note the distance between the upper part of the great trochanter and the distal tip of the prosthesis to be implanted. Add 1 cm to this distance.
- Select the cement restrictor matching the sizing trial.
- Unlock the sizing trial from the inserter. Attach the restrictor to the inserter.
- Introduce the cement restrictor into the medullar cavity to the previously predetermined distance.
- Note: a distance of 1cm is advised between the distal tip of the prosthesis and the restrictor.
- Unlock the set screw.
- Pull the handle upwards to leave the restrictor in situ. The inserter may now be withdrawn. The cement restrictor is now positioned within the canal at the predetermined distance.
- Introduce the acrylic bone cement.

Caution : the instrument set must not be forced.

WARNINGS AND PRECAUTIONS FOR USE
 Read instructions for use carefully prior to use of **CEMSTOP®** device. This device contains a biological substance from animal origin (gelatin from pigskin). Safety of this substance has been verified, no additional caution or precautions for use are required unless patient is allergic. Ignoring the instructions for use may lead to potential undesirable effects. Safety of the **CEMSTOP®** degradation products has been verified, no additional caution or precautions for use are required unless patient is allergic. Ignoring the instructions for use may lead to potential undesirable effects.

For an optimum use of **CEMSTOP®**, it is essential to perform a thorough pre-operative review of the patient before the intervention in order to confirm the indication and plan the surgical technique.

Conservation
- It is strictly forbidden to re-sterilise the product. This product is supplied sterile unless package has been opened or damaged.
 - This device is packaged and sterilised for single use only. Do not reuse, reprocess or resterilise. Reuse, reprocessing, or resterilisation may compromise the structural integrity of the device and/or lead to device failure which may result in patient injury or illness. Also, reprocessing or resterilisation of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another.
 - Before use, check the protective wrapping carefully to ensure that it has not been opened nor damaged in a way that could affect its sterility.
 - Do not use this product after the expiration date printed on the package.

Use
 - Examine the product visually to identify any defects such as crack or deformation. Do not implant products with defects.
 - When removing the product from its wrapping, be sure to follow asepsis rules.
 - If packaging is unintentionally opened before use or damaged, do not use the product.
 - Special attention should be given to the choice of the size of the cement restrictor to avoid risks of cement flowing down or tearing of the cement restrictor.

INFORMATION TO THE PATIENT
 The patient should be informed by the operator of potential consequences of factors mentioned in the paragraphs « contra-indications » and « side effects, and adverse events », meaning those liable to hinder the success of the operation, as well as the possible complications which may arise. The patient should also be informed that this device contains a biological substance of animal origin (gelatin from pigskins). He should also be informed of measures to be taken to diminish possible consequences of these factors.

The patient must be informed by the operator that after implantation, this device does not require any precautions or specific measures regarding the external environment conditions (eg : security portals, cooking induction plate) or other medical exams (eg : MRI).

STORAGE
 Each product must be stored unopened in its original packaging in a dry place and at a temperature below 25°C.

RECOMMENDATIONS FOR DISPOSAL
 Comply with local regulations in force related to handling and disposal of the medical device and its packaging.

IMPLANT CARD
 As part of traceability of its devices, the manufacturer is providing

for the patient an implant card on which a label specific to the device must be applied. The implant card as well as the label are located in the device box. Patients shall keep this implant card at all times in their medical records.

Instruction for completion:
 •  ?: Patient identification. To be filled by the healthcare institution.
 •  ?: Date of implantation. To be filled by the healthcare institution.
 •  ?: Name and address of the healthcare institution. To be filled by the healthcare institution.
 • "IMPLANT CARD" sticker with device information. To be affixed by the healthcare institution.

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